EndoArm

Product Manual

Fillauer®

Table of Contents

English	3
Spanish.....	16
French	28
German	40
Italian	52
Swedish	64
Norwegian.....	76
Finnish	88
Danish	100
Dutch	112

Intended Use/Purpose

The EndoArm system is a revolutionary arm and elbow, modular solution for passive transradial and transhumeral prostheses. The EndoArm Flexion/Friction wrist offers adjustable friction, flexion and rotation to allow the user to easily re-position the hand for optimization of his/ her task.

The EndoArm Elbow offers the user infinite locking positions over its 150° ROM along with an adjustable lift assist. The EndoArm components offer a lightweight modular system to provide the user with the functions needed for everyday living.

The connecting threads come standard in M12*1,5mm. Some are available with a ½"-20 thread as well.

EndoArm Performance Characteristics

- Adjustable function. Add features as needed
- Adapted for 20mm aluminum tubing
- Connectors are adapted for ½"-20 and M12x 1,5 threads internal and external
- Length: adaptable depending on what parts are combined
- Material: Aluminum, brass, POM and steel

EndoArm Elbow Performance Characteristics

- Shortening of the forearm when flexing the elbow to reduce the wear on the cosmetic foam covering the finished prosthetic
- The elbow unit comes with the lock slide set so the elbow can be locked by pressing the button from the lateral side.
- This can be changed in the workshop if needed.

- The proximal attachment features a M12*1.5mm internal thread or using a M8 bolt attachment to the socket (M8 bolt not included).
- Tighten or loosen the brass lift assist adjustment nut to increase or decrease the lift assist as needed.
- Complete modularity with our other EndoArm Products
- Length: 5.5 in. (140 mm)
- Weight: 4.6 oz. (132 g)
- Load: 8 kg at hand

Storage and Handling

It is recommended that EndoArm System be stored in a cool, clean, dry environment away from harsh chemicals (chlorine, acids, acetone, etc.).

Keep the EndoArm components clean and dry, lubricate moving parts with a thin oil when needed.

Warnings and Precautions

WARNING: This device has been designed and tested for activities of daily living and not for abnormal activities such as extreme sports. To avoid risk of injury, handle your prosthesis and its components with care. If your prosthetic components are subjected to abnormal stress (such as a fall) contact your prosthetist/orthotist to have the components inspected for any damage.

NOTICE: The device is intended for single patient use only.

NOTICE: The distal attachment for the EndoArm Elbow is a 20mm tube clamp, standard to our EndoArm Tube. Please note that this tube cannot be shorter than 65mm due to the lift assist.

NOTICE: Attachment, adjustment, alignment, and delivery of this device must be performed by or under the direct supervision of a qualified prosthetist/orthotist. Unless stated in this manual, any such activities should not be attempted by the user and will potentially void the device warranty.

NOTICE: Abnormal or improper environmental conditions will lead to malfunction and damage of the prosthesis and is not covered under the warranty of the device. This prosthetic component must not be subjected to dust/debris, harsh chemicals, abrasives, vibration, activities which would damage the biological limb, or prolonged exposure to extreme temperatures ($<-15^{\circ}\text{C}$ or $>50^{\circ}\text{C}$) ($<5^{\circ}\text{F}$ or $>122^{\circ}\text{F}$). Do not allow debris or liquids to remain in the prosthesis/orthosis and its components during use.

NOTICE: The ability of an upper limb prosthetic device user to drive should be determined on a case-by-case basis by a specialist. Contact your local governing authorities regarding any driving restrictions or limitations. Myoelectric devices should be turned off and/or have electric movement disabled and body-powered devices should not rely on cable tension for grasp control if the user has been cleared to drive with the prosthesis.

NOTICE: When creating a cosmetic cover for the EndoArm Elbow, ensure that no material can get in between the lock/unlock slide button and the elbow housing as this will affect the elbows performance and reliability.

Specifications and Preparations Before Use

EndoArm Tube

631028330 – EndoArm Tube 0.8 x 5.9 in. (20 x 150 mm)

Cut the tube to the required length with either a pipe cutter for aluminum or a hacksaw. Ensure that the tube keeps its circular cross-section. Deburr edges before use.

Instructions for EndoArm-Friction Flexion Rotation Wrist

631028300 - Flexion / friction / rotation wrist

Assemble unit with a small amount of Loctite 242 or similar, applied to the M5X20 thread to ensure the chosen friction setting of flexion/ extension resistance.

Tighten the M3X16 bolt to ensure that the tube is secured in place. Do not exceed 2Nm torque.

The M12X1,5 mm threaded stud is secured preferably to the chosen hand with a suitable adhesive.

EndoArm - Tube Clamp Adapter

631028310 - Tube Clamp, external thread 1/2"-20

631028311 - Tube Clamp, internal thread 1/2"-20

631028340 - Tube Clamp, external thread M12 x 1.5

631028341 - Tube Clamp, internal thread M12 x 1.5

Tighten the M3X16 bolt to ensure that the tube is secured in place. Do not exceed 2Nm torque, as this might damage the bolt. The M12X1,5mm or 1/2-20" threaded stud is secured preferably to the chosen hand with a suitable adhesive.

EndoArm - Elbow

631028349- EndoArm Elbow, Left

631028350- EndoArm Elbow, Right

631028351 -Lock Plate x 2 and Lock Rod (spare parts)

631028353- TH Friction Adapter (Addon)

Tighten the bolt M4x18 bolt at the tube clamp to approximately 2Nm and apply a suitable thread lock before final delivery.


Change lock direction EndoArm Elbow

The Slidelock is set by default so that the locking motion occurs from the lateral side of arm.

This is so that the user can use its surroundings to assist in locking the elbow


Step 1

Remove the lift assist adjustment nut (1) and lift assist spring (2). Then remove screws (3) and (4) so that the lid (5) can be removed.


Step 2

Unscrew the Lock rod (6) and pull it out so that it goes clear of the Slide lock assembly (7).

Step 3


Lift the Slide lock assembly (7) straight up. Try to leave the Lock plates and Lock plate springs (8) in place.

Step 4


Rotate the Slide lock assembly (7) perpendicular to its longitudinal axis so that the buttons profile still matches the out- side contour of the Elbow tube clamp (9).

Insert the Slide lock assembly (7) into the Elbow tube clamp (9) and reinsert the Lock plates and Lock plate springs (8)

Step 5

Reinsert the Lock rod (6) through the unit and reattach to the brass Lock rod seat (10). Apply a small amount of Loctite 243, or similar, to the thread if needed.

Step 6

Reattach the parts removed in Step 1. Test the functionality and motion of the elbow before reattaching to the prosthesis.


Consumable Components

Replace Lock plate and Lock rod EndoArm Elbow

631028351 -Lock Plate x 4 and Lock Rod (spare parts)


Step 1

Remove the lift assist adjustment nut (1) and lift assist spring (2). Then remove screws (3) and (4) so that the lid (5) can be removed

Step 2


Unscrew the Lock rod (6) and pull it out so that it goes clear of the Slide lock assembly (7) lock plates and springs (8).

Step 3


Use a pair of tweezers to remove the lock plates (8) and springs (9). Insert the parts from Replacement kit 63128351 instead.

A pair of tweezers will help you reinsert the lock plates. Use the same tweezers to insert the compressed springs into the elbow tube clamp housing (10).

Step 4

Reinsert the Lock rod (6) through the unit and reattach to the brass Lock rod seat (10). Apply a small amount of Loctite 243, or similar, to the thread if needed

**Step 5**


Reattach the parts removed in Step 1. Test the functionality and motion of the elbow before reattaching to the prosthesis. Apply some low viscosity lubricant to the lock rod if needed.

Attaching TH Friction Adapter EndoArm Elbow

631028353- EndoArm Elbow TH Friction Adapter (Addon)

Step 1

Attach the Nylon Washer 25x2,5x10,5 (3) to the stud of the EndoArm Elbow TH Friction Adapter (1).


Step 2

Insert into the Endoarm Elbow


Step 3

Attach M6x40 Dome Bolt (5) with Nylon Washer 12x1,6x6,4 (4) through the EndoArm Elbow and thread into the friction adapter

Step 4

From the top attach the M6 Lock nut (6), but do not tighten yet.

Step 5

Using a 4mm Hex key tighten the M6 Dome Bolt to give the required friction for the transhumeral assembly.

Lock the friction setting by tightening the M6 lock nut using a 10mm wrench while using a 4mm hex key to hold the bolt in place.

Compatibility (Combination of Devices)

Fillauer has tested (ISO 22523:2006) and recommends the use of standard, adult, endoskeletal components from Fillauer with all passive transradial and transhumeral prostheses. Components from other manufacturers may or may not be compatible. Failure due to use of other manufacturers' products is not covered under warranty.

The following components are intended to be used in conjunction with the EndoArm System:

- 631028349 - EndoArm Elbow, Left
- 631028350 - EndoArm Elbow, Right
- 631028401 - EndoArm Kit with Flexion/ Friction/ Rotation Wrist
- 631028402 - EndoArm Kit
- 631300100 - Cosmetic Foam for TR
- 631028300 - Flexion / friction / rotation wrist
- 631028306 - Center screw, M5 x 20 mm
- 631028310 - Tube Clamp, external thread 1/2"-20
- 631028311 - Tube Clamp, internal thread 1/2"-20
- 631028320 - Lamination adapter, M12 x 1.5 x 20 mm
- 631028330 - Tube 0.8 x 5.9 in. (20 x 150 mm)
- 631028340 - Tube Clamp, external thread M12 x 1.5
- 631028341 - Tube Clamp, internal thread M12 x 1.5
- 631028342 - Clamp screw, M3 x 16 mm

Disposal/Waste Handling

The product must be disposed of in accordance with applicable local laws and regulations. If the product has been exposed to bacteria or other infectious agents, it must be disposed of in accordance with applicable laws and regulations for the handling of contaminated material.

User Information

The providing health care professional must review the following information directly with the user:

Attachment, alignment, and delivery of this orthotic/prosthetic device must be performed by or under the direct supervision of a qualified prosthetist. Any adjustment or modifications should be done by the clinician and not by the user.

Warnings and Precautions For The User

WARNING: This device has been designed and tested for activities of daily living and not for abnormal activities such as extreme sports. To avoid risk of injury, handle your prosthesis and its components with care. If your prosthetic components are subjected to abnormal stress (such as a fall) contact your prosthetist/orthotist to have the components inspected for any damage.

NOTICE: Abnormal or improper environmental conditions will lead to malfunction and damage of the prosthesis and is not covered under the warranty of the device. This prosthetic component must not be subjected to dust/debris, harsh chemicals, abrasives, vibration, activities which would damage the biological limb, or prolonged exposure to extreme temperatures (<-15C or >50C) (<5F or >122F). Do not allow debris or liquids to remain in the prosthesis/orthosis and its components during use.

NOTICE: The ability of an upper limb prosthetic device user to drive should be determined on a case-by-case basis by a specialist. Contact your local governing authorities regarding any driving restrictions or limitations. Myoelectric devices should be turned off and/or have electric movement disabled and body-powered devices should not rely on cable tension for grasp control if the user has been cleared to drive with the prosthesis.

Consultation of Health Care Professional

The user should monitor their prosthesis daily and contact their health care professional if they experience changes in performance of the device or hear irregular noise from it.

Care and Maintenance

Notice: Handle your prosthesis and its components with care and keep them free from debris and moisture

Serious Incidents

In the unlikely event a serious incident occurs in relation to the use of the device, users should seek immediate medical help and contact their prosthetist, local competent authority and Fillauer at the earliest possible convenience. Clinicians should at any time contact their local Fillauer representative and local competent authority immediately in the event of any device failure.

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the Instructions for Use. The device must only be operated with tested modular components. The manufacturer is not liable for damage caused by component combinations that were not authorised by the manufacturer.

Warranty

This product has a 12 month warranty against manufacturer defects.

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© 2023 Fillauer LLC