Comfil Material

Product Manual

Fillauer®

Table of Contents

English	3
Spanish.....	11
French	19
German	27
Italian	35
Swedish	43
Norwegian.....	51
Finnish	59
Danish	67
Dutch	75

Intended Use/Purpose

The COMFIL material combines the high localized strength of carbon composite with the modeling properties of thermoplastics. Available in unidirectional and bidirectional, the carbon fibers are interwoven with a revolutionary thermoplastic fiber that allows the material to be heated and formed at lower temperatures. This thermoplastic material also acts as a high-strength matrix for active performance and flexibility. It is available as sheets, insoles and bars.

Sheets

Thermo Formable Composite UD Form

PRODUCT NUMBER	DESCRIPTION	GODD FOR:
700 010 101	UD Form, Carbonfiber 1000 x 600 x 1,6 mm	Good for Bars and other rigid parts
700 010 103	UD Form, Carbon & Glassfiber 1000 x 600 x 2,55 mm	Good for Bars and rigid insoles
700 010 126	UD Form, Carbonfiber 1000 x 600 x 0,7mm	Good for Reinforcement

Thermo Formable Composite BI Form

PRODUCT NUMBER	DESCRIPTION	GODD FOR:
700 010 121	BI Form, Glassfiber black 1000 x 600 x 0,75 mm	Opening parts
700 010 122	BI Form, Glassfiber black 1000 x 600 x 1,5 mm	Flexible toe-off parts
700 010 123	BI Form, Carbon & Glassfiber 1000 x 600 x 1,5 mm	Rigid foot plate
700 010 124	BI Form, Carbonfiber 1000 x 600 x 0,7 mm	Difficult shapes
700 010 127	BI Form, Carbonfiber 1000 x 600 x 0,5 mm	Difficult shapes
700 010 128	BI Form, Glassfiber 1000 x 600 x 0,7 mm	Difficult shapes
700 010 130	BI Form, Aramid 1000 x 600 x 0,9 mm	Flexible parts

Insoles, AFO Kit & Bars

Thermo Formable Composite Insole

PRODUCT NUMBER	DESCRIPTION	PAIR
700 030 121	Insole Flexibel, Glassfiber Black 1,5 mm	Size 35 - 37
700 030 122	Insole Flexibel, Glassfiber Black 1,5 mm	Size 38 - 40
700 030 123	Insole Flexibel, Glassfiber Black 1,5 mm	Size 41 - 43
700 030 124	Insole Flexibel, Glassfiber Black 1,5 mm	Size 44 - 46
700 030 141	Insole Rigid, Carbon & Glassfiber 2,5 mm	Size 35 - 37
700 030 142	Insole Rigid, Carbon & Glassfiber 2,5 mm	Size 38 - 40
700 030 143	Insole Rigid, Carbon & Glassfiber 2,5 mm	Size 41 - 43
700 030 144	Insole Rigid, Carbon & Glassfiber 2,5 mm	Size 44 - 46

AFO Kits

PRODUCT NUMBER	DESCRIPTION	LENGTH
700 020 111	Small	A - 148 mm B - 150 mm C - 600 mm
700 020 112	Medium	A - 168 mm B - 185 mm C - 600 mm
700 020 113	Large	A - 195 mm B - 215 mm C - 1000 mm

Bars

PRODUCT NUMBER	DESCRIPTION
700 020 321	UD Form, Carbon-& Glassfiber 600 x 20 x 2,55 mm
700 020 331	UD Form, Carbonfiber 600 x 20 x 1,6 mm
700 040 420	UD Form, Carbonfiber 600 x 20 x 3,2 mm

Clinical Benefits

COMFIL gives you the possibility to produce stronger and thinner orthotics and prosthetics compared to non-reinforced thermoplastics. The thermoplastic matrix makes it possible to make adjustments to your device after fabrication. After fabrication adjustments are not possible with a thermoset matrix or prepreg material.

Performance Characteristics

The COMFIL material can be simply custom formed to your mold with a heat gun or by placing the material into an oven to achieve molding and extra reinforcement materials can be simply heated and applied to customize your brace for the patient's needs. Your custom lay-up is then placed into a standard oven under vacuum at a maximum of 230 degrees Celsius for 30- 45 minutes to consolidate the lay-up into one piece for an individual orthotics or prosthetic solution. The COMFIL material is also fully heat remoldable and allows adjustable fitting just like thermoplastics and unlike prepreps. The material comes in unidirectional and bidirectional reinforcement for different structural properties with three different thicknesses and can be cut simply with scissors.

Material

- Thermoplastic
- Carbon fiber
- Glass fiber

Storage and Handling

Do not place this product in direct sunlight for long periods of time. When not in use, store in a cool, clean, dry, and dark place.

Warnings and Precautions

WARNING: This product contains continuous strands of Carbon fiber.
Possible hazards of carbon fiber are:

- Mechanical irritation (itching)
- The formation of respirable dust and non-respirable filaments
- Extremely rare possibilities of allergy.

Exposure controls and personal protection

NOTICE: It is recommended your work area is adequately ventilated, especially in confined areas.

Required personal protection equipment:

- Respiratory protection: Breathing apparatus in the event of high concentrations.
Short term: filter apparatus, Filter P1.
- Hand protection: Butyl rubber, > 120 min (EN374)
- Eye protection: Safety glasses
- Skin and body protection: lightweight protective clothing
- Hygiene measures: General industrial hygiene practice. Regular cleaning of equipment, work area and clothing.

Specifications and Preparations Before Use

Preparation


Plaster Preparation:

- Dry the plaster
- Cover the plaster with a cotton stockinette
- Draw out the orthoses design

Material

- Cut the parts of the various Comfil materials you want to use. Make sure you get the right strength and flexibility through fiber direction
- The UD material should be cut in the length direction of the fibers.
- The BI material should be cut at a 45° angle on the sheet to allow for easy shaping and flexibility.
- Cut the pieces as close to the finished models you want. Anything to reduce grinding after completion
- Grind the material with long, sweeping movements (avoid overheating the material when grinding)

Manufacturing

Heat the Comfil material with a heat gun at the highest temperature in a sweeping motion until the material becomes ductile. Shape the material into the basic shape with light pressure. Finish shaping by pressing the material firm against the plaster.

Then proceed with the next part. Overlapping material for maximal sufficient adhesion and strength. Begin to shape the material by affixing it on the already formed part. Make sure there is no air between the different elements and between the material and plaster.

Continue section by section until the orthosis lay up is finished. Before the orthosis is laminated in the oven grind the edges smooth. After grinding, if necessary, try the orthosis on the patient. Remember that the orthosis is not yet at full strength.

When all parts are in place, take it off the plaster. Dress the plaster with a release film. Use heat-resistant tape to attach the release film. Place the orthosis back on the plaster, be careful not to damage the film.

Release Film - 700 000 023

Tape - 700 000 015

To obtain a shiny surface on the orthosis stretch a silicon sheet over the orthosis, stretch the material enough to prevent wrinkles. Use heat-resistant tape to attach the material. Then place a cotton stockinette over the silicone, as a breather for the vacuum.

Silicone sheet - 700 000 016

Tape - 700 000 015

Cotton stockinette - 636 508 020

Prepare the vacuum bag by sealing one end with sealing tape. Prepare the vacuum stand by wrapping sealing tape under the vacuum hole. Cover the plaster with the vacuum bag and then seal it by tying the bag against the sealing tape. Pull vacuum on the bag. For best results you need .08 bars or higher vacuum. Place the plaster in the oven for approx. 40 min at 230o Take out the plaster and let it cool under vacuum. Important: the vacuum must be on as long as the material is hot. The plaster can be cooled rapidly by placing it in cold water.

Sealing tape - 700 000 014

Vacuum Bag - 700 000 012

Requirements Of Competencies

NOTICE: Attachment, adjustment, alignment, and delivery of this device must be performed by or under the direct supervision of a qualified prosthetist/orthotist. Unless stated in this manual, any such activities should not be attempted by the user and will potentially void the device warranty.

Combinations Of Materials

Comfil Material can be used as the whole material but also as a reinforcement in orthotics and prosthetics.

Disposal/Waste Handling

The product must be disposed of in accordance with applicable local laws and regulations. If the product has been exposed to bacteria or other infectious agents, it must be disposed of in accordance with applicable laws and regulations for the handling of contaminated material.

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© 2023 Fillauer LLC