
A technical line drawing of a jaw crusher mechanism, showing the main frame, the movable jaw, the fixed jaw, and the toggle plate assembly. The drawing is rendered in white lines on a dark blue background.

Fillauer TRS, Jaws (Jaws, Jaws XFS)

Product Manual

Fillauer®

Table of Contents

English	3
Spanish.....	11
French	19
German	27
Italian	36
Swedish	43
Norwegian.....	51
Finnish	59
Danish	67
Dutch	75

Intended Use/Purpose

JAWS and JAWS XFS (Extra Force Spring) are heavy duty, voluntary opening (VO), body powered, terminal devices, which operate universally as either right or left-handed. JAWS are waterproof and are truly a “cross-over” type of prosthetic device capable of enhancing control over tools, handles and equipment as well as very functional for vehicle control of work and recreational vehicles such as forklifts, tractors, ATVs, snowmobiles, motorbikes and cycles and jet ski type watercraft. JAWS are a perfect, functional compliment to users of VO split hook technology providing a range of performance and control not possible with that type of technology. JAWS mount easily to any USA standard prosthetic wrist.

***JAWS XFS is only intended for users who need extremely strong gripping forces. JAWS XFS should be used with extreme caution, especially in the highest force position, because “releasing the grip” has an extremely high threshold. Additionally, if operating the JAWS XFS with a traditional cable system it may cause premature wear and possible failure of the cable system, due to the extra force required to operate the TD.**

The device is intended for single user/patient use only.

Performance Characteristics

JAWS can be operated in three ways. Cable operation, Manual “pull open” using a built-in, “trigger” feature or Push-On/Pull-Off (Snap ON-OFF) feature. JAWS models have no lock system for safety. JAWS have a manual, shift mechanism, providing 4 gripping force levels. JAWS are equipped with large compliant but tough synthetic rubber gripping surfaces for optimum control over handles and cylindrical objects. JAWS have a manually adjusted, friction, wrist pivot system that provides the user with greater range of motion control over objects reducing elbow and shoulder joint movement and compensation.

Indications/Contraindications

General:

Upper Extremity Terminal Devices-End Effectors:

Prosthetic devices are tools designed to provide or replicate certain aspects and functions of the human hand including anatomical realism. All these devices have benefits, limitations, ranges of function and liabilities that need to be understood by the user-consumer.

The basis for any reliable functional outcome in using these devices is premised upon having a well-designed, reliably suspended, comfortable, functional prosthetic limb, which optimizes the user's remaining physical capabilities, including range of limb motion and strength. The remaining limb itself, is a very basic limiter of function and performance. Typically, the longer the limb the more functional capacity over a prosthesis will be achieved, assuming the limb's muscular and skeletal framework are not compromised beyond simply hand absence. Loss of muscle tissue, permanent nerve damage and phantom limb sensations-pain can all impact the user's ability to perform with a prosthesis. In general, a user with a limb absence where less than 25% of the humerus remains, will **not** be a good candidate for using upper extremity prosthetic technology successfully. Additionally, the user's cognitive acuity and capability are important in understanding the ranges of function and the specific limits of function, while controlling such prosthetic tools, to optimize their performance and avoid injury.

Finally, all prosthetic terminal devices, have inherent dangers of entanglement or engagement where release can be compromised because of their physical exterior design, unique unforgiving materials, and inanimate lack of "feel". Wearing an upper extremity prosthesis does involve risk! *Training and therapy are always recommended when using a prosthesis, especially when using new or unique technology or changing prosthetic system operative controls.*

Specific:

Body Powered Devices: Prehensor TDs.

Body powered prosthetic devices using a shoulder harness and cable for operation require a certain level of physical coordination and strength. Additionally, use of such devices requires a sound cognitive basis that provides the user with the understanding of the values of the technology, how it operates and its limitations and the potential for both performance and self-injury. Body powered prosthetic devices require practice. More practice will provide better functional outcomes in performance. Certain devices have capabilities beyond that of the human hand regarding their ability to withstand impact, withstand load and resist environmental elements such as severe heat, fire, freezing cold temperatures and caustic chemicals. Other devices are inferior to the human hand regarding such exposures. The user must understand and take advantage of the aspects of each type of technology they choose to use.

Storage and Handling

It is recommended that JAWS (or prosthetic/orthotic components) are stored in a cool, clean, dry environment away from harsh chemicals (chlorine, acids, acetone, etc.).

Warnings and Precautions

NOTICE: The ability of an upper limb prosthetic device user to drive should be determined on a case-by-case basis by a specialist. Contact your local governing authorities regarding any driving restrictions or limitations. Voluntary Opening body-powered devices can rely on cable control or external spring/elastic band tension for grasp control, even if the user has been cleared to drive with the prosthesis.

CAUTION: Abnormal or improper environmental conditions will lead to malfunctioning and damage of the prosthesis and is not covered under the warranty of the device. This prosthetic component must not be subjected to dust/debris, liquids other than fresh water, abrasives, vibration, or activities which would damage the biological limb. Do not allow debris or liquids to remain in the prosthesis and its components during use. Rinse the wrist with fresh water and dry immediately after exposure.

CAUTION: The JAWS is waterproof to 1 meter; however, if JAWS is submerged, it should be rinsed with fresh water and **dried** immediately to remove salt, chlorine, or debris.

Qualified provider

Attachment, adjustment, alignment, and delivery of this device must be performed by or under the direct supervision of a qualified prosthetist. Unless stated in this manual, any such activities should not be attempted by the user and will potentially void the device warranty.

Specifications and Preparations Before Use

Product Name: Jaws Voluntary Opening Prehensor

Specifications

Length	5.0 in. (12.7 cm.)
Width	2.1 in (5.3 cm.)
Height	3.9 in. (9.9 cm.)
Weight	15.5 oz. (440 gm.)
Color	Aluminum & Black

Applications

Age	Older teenagers and adults
Limb Description:	All models: Trans-radial & mid-length - long Trans-humeral

Preparation

Consider the following factors as you determine the most appropriate technology and prosthesis design for your client:

1. Does the user understand the function and safety features of the device?
2. What other terminal devices will be used with this prosthesis?
3. In what types of environments will the device(s) be used?
4. In what positions will the person use the device?
5. How will the user switch terminal devices? Will quick-disconnect be helpful?
6. For what other activities might this device be helpful?
7. Will the person benefit from additional training in the use of the device?

Installation

JAWS may be installed in any Fillauer TRS Omega Wrist, and any Fillauer LLC wrist unit with a ½-20 thread. Follow the instructions provided with the wrist unit for best results.

Cabling

The “thumb” has been designed with a receiver for a 9/32” terminal ball and a 3/16” terminal ball. A Binary Swivel, or triple swivel with either a 3/16” or 9/32” terminal ball should be selected with a cable connection that matches the cable used in the prosthesis. Cable routing should assure a direct line of pull that minimizes bends in the cable which could result in excess cable friction or failure.

Tension of the JAWS gripping force is adjusted by a manual shift mechanism, providing four gripping force levels.

Compatibility

Important CLINICAL NOTE & FEATURE about all body powered VO and VC terminal devices

Body powered terminal devices develop control through the biomechanical manipulation of body forces that are captured through a shoulder harness system and transmitted through a cable to the terminal device. The function, operation, ease of use and optimization of body power in a cable powered terminal device can be improved by the addition of an “In-Line” cable modulation and cable locking system, such as the Fillauer TRS SURE-LOK technology. Such technology allows the prosthetic user to grasp objects then reduce or totally relieve tension on the cable system. This function increases the range of object control motion available to the user and reduces fatigue, while combatting the development of “over-use” syndrome.

JAWS have been tested with and are recommended for use with Fillauer TRS Omega Wrists, and Fillauer LLC wrists that have a ½-20 internal thread. They may be used with any equivalent ½-20 threaded wrist units; however, damage caused by other manufacturers wrist units is not covered under warranty of this device.

Care and Maintenance

General:

Upper Extremity Terminal Devices-End Effectors:

Prosthetic devices should be used and treated as quality tools, regularly inspected for damage, deterioration and or excessive wear. The devices are waterproof and or highly water resistant and constructed of materials that will not rust or corrode or are resistant to such exposure. They can be immersed in water and washed. If truly mechanical they should be dried, as thoroughly as possible and light lubricants can be applied to moving mechanical linkages as needed. Soap and water cleansing or rinsing are recommended, especially after operation in environments with a lot of dirt, dust, grit, oils or solvents. If used in ocean water, cleaning as soon as possible in fresh water is highly recommended with all devices. Long term exposure in intense direct sunlight or heat outside of normal use will damage prosthetic devices that incorporate synthetic rubber in their construction or coverings.

Specific:

Body Powered Devices: Prehensor TDs.

These devices all employ certain mechanics for their operation. Regular daily cleaning-washing is reasonable and recommended. Treat these devices as if they were your hand from a sanitation perspective. Use a light lubricant, as might be required, on moving parts and linkages. Do not use heavier viscous oils or lubricants. Dry film lubricants with silicone are fine as well. Regularly inspect these types of devices for wear, loose fasteners or fittings or other aspects of deterioration. If damage or wear is observed, then have a qualified prosthetic technician complete the repair, replace parts, etc. to ensure that the device performs optimally and safely.

Disposal/Waste Handling

The product must be disposed of in accordance with applicable local laws and regulations. If the product has been exposed to bacteria or other infectious agents, it must be disposed of in accordance with applicable laws and regulations for the handling of contaminated material.

All metal components may be removed and recycled at the appropriate recycling facility.

User Information

The providing health care professional must review the following directly with the user:

Preparation for User

Consider the following factors as you prepare for this device:

1. Are you able to don and to doff the device independently?
2. Are you able to attach/detach the cable from the device?
3. Are you able to switch terminal devices independently?
4. Are you able to access services from a skilled occupational therapist for additional training if needed?

Warnings and Precautions for the user

NOTICE: The ability of an upper limb prosthetic device user to drive should be determined on a case-by-case basis by a specialist. Contact your local governing authorities regarding any driving restrictions or limitations. Voluntary Opening body-powered devices can rely on cable control or external spring/elastic band tension for grasp control, even if the user has been cleared to drive with the prosthesis.

CAUTION: Abnormal or improper environmental conditions will lead to malfunctioning and damage of the prosthesis and is not covered under the warranty of the device. This prosthetic component must not be subjected to dust/debris, liquids other than fresh water, abrasives, vibration, or activities which would damage the biological limb. Do not allow debris or liquids to remain in the prosthesis and its components during use. Rinse the wrist with fresh water and dry immediately after exposure.

CAUTION: The JAWS is waterproof to 1 meter; however, if JAWS is submerged, it should be rinsed with fresh water and **dried** immediately to remove salt, chlorine, or debris.

Serious Incidents

In the unlikely event of a failure resulting in a fall and/or injury, seek immediate medical help and contact your prosthetist at the earliest possible convenience.

Warranty

This product has a 12-month warranty against manufacturer defects

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© 2021 Motion Control, Inc. 1910013 Rev H 06-29-2021
Motion Control, Inc.