Mechanical Lock Gripping Devices

Product Manual

Fillauer®

Table of Contents

English	3
Spanish.....	14
French	26
German	38
Italian	50
Swedish	62
Norwegian.....	74
Finnish	86
Danish	98
Dutch	109

Intended Use/Purpose

All the products listed above have the same basic function in that they are manually engaged and disengaged, mechanical, locking, gripping terminal devices. They can be operated as either right or left models. They are not cable operated and, in the USA, are most often classified as “Activity Specific” terminal devices. They are designed to capture and hold objects with a clamp, flexible strap or ratchet closure strap system for optimal use with cylindrical and semi-cylindrical shaped objects and handles. They all provide stable, strong and secure control over the objects engaged but all require a manual “release” of a mechanism or strap to disengage from those objects. All these devices mount and connect easily to any USA standard prosthetic wrist.

These devices are intended for single user/patient use only.

Performance Characteristics

BLACK IRON LITE, BLACK IRON LITE PIVOT, BLACK IRON MASTER, BLACK IRON TRAINER, BLACK IRON TRAINER PIVOT

All these terminal devices were primarily designed for resistance exercise using barbells, dumbbells (Free weights), as well as other types of exercise equipment where handles are used and engaged to perform and complete the activity. The PIVOT models are slightly longer but provide an additional range of motion to help mitigate or eliminate excessive compensatory movements at the elbow and shoulder for more natural biomechanics throughout the exercise. The BLACK IRON LITE (BIL) and BLACK IRON LITE PIVOT are lighter-duty devices designed primarily for aerobic dumbbell exercise.

They are not recommended for weight resistance activities in excess of 10 pounds (4.5 KG). They incorporate a ratchet strap system to secure the weights into the terminal device. The BLACK IRON TRAINER (BIT) AND BLACK IRON TRAINER PIVOT models, which weigh about 1 pound (440 gm), are the most popular designs capable of handling loads in excess of 400 pounds (~200 kg). These models will meet the needs of almost all persons with a hand absence wanting to weight train. The BIT and BIM both use a swing up “gate” closure and large bolt and wing nut assembly to capture handles and securely lock in the weights to the prosthesis. The BLACK IRON MASTER (BIM), which weighs 25 ounces (~700 gm.), is a professional level weight training device capable of controlling loads of over 1600 pounds (~700kg). It is recommended for use only by professional body builders and other athletes who seek exceptionally high performance and load capacity in their exercise training regime.

FLEXIBLE THING

The Flexible Thing is a terminal device designed primarily for water sports like water skiing, where bilateral control over a “tow handle” is required but where safety necessitates the ability to release the handle instantly during a mishap or fall to avoid injury. The device requires the “over-grasping” of a strap, which is wrapped around the handle, by a “sound” hand, to keep it engaged onto any tow handle. The Flexible Thing is also equipped with a special buckle into which the strap can be woven for a multitude of other tasks where control over cylindrical handles typical in shovels, rakes, etc. is required for yard work or carrying heavy handled objects. This feature provides the Flexible Thing to be categorized as a true “cross-over” type of terminal device because of its multi-use applications and capabilities.

HAMMERHEAD, HAMMERHEAD PIVOT

These terminal devices were primarily designed to excel in performance in the activity of kayaking but have been applied to canoeing as well with more limited success. The Hammerhead uses a solid polymer body and strong elastic strap which wraps over and secures the kayak paddle shaft into the device. A cylindrical handle in the end of the strap allows the strap to be manually snapped down and locked into position over a pedestal. The same handle provides for the ability to manually grasp and release the strap thereby freeing the paddle from the control of the device. The Hammerhead operates either right or left and replicates the biomechanical motions of the hand, wrist and forearm used to propel a kayak or similar water craft. The Hammerhead Pivot model is longer and adds additional radial-ulnar movement at the wrist that enhances performance and helps mitigate or eliminate excess compensatory movements at the elbow and shoulders while paddling. Both devices are totally waterproof and can be used in both fresh and salt water environments.

ISHI, ISHI HD, MULTI-D, MULTI-D PIVOT

These products are truly “cross-over” terminal devices and can be used for a wide variety of holding tasks and are especially created for holding cylindrical shaped handles and objects. They operate as either right or left models. Both the ISHI and Multi-D utilize flexible grasping “fingers” or “mandibles” biased together with high strength, quick release, ratchet strap, closure mechanisms. The ISHI and ISHI HD were originally conceived and designed for holding a bow or fishing rod but they function very well in holding other tools and implements in a narrow diameter range up to ~2 inches (~5 cm). The MULTI-D is a larger version of the ISHI design but with broader mandible gripping surfaces and a heavier capacity ratchet strap closure system for heavier duty use and activities. The Multi-D can handle objects of bigger diameter out to ~3 inches (7.6 cm). The Multi-D Pivot is longer, slightly heavier and provides radial-ulnar or flexion-extension types of wrist movement depending on how it is mounted on the prosthetic wrist. Both these devices are waterproof. The Multi-D has been used for vocational tasks like using drills and hammers but can also “crossover” and be used for holding kayak, canoe and standup style, paddle board paddles and for controlling boat oars, which are large in diameter than most paddle-style handles and grips.

Indications/Contraindications

General:

Upper Extremity Terminal Devices-End Effectors:

Prosthetic devices are tools designed to provide or replicate certain aspects and functions of the human hand including anatomical realism. All of these devices have benefits, limitations, ranges of function and liabilities that need to be understood by the user-consumer. The basis for any reliable functional outcome in using these devices is premised upon having a well-designed, reliably suspended, comfortable, functional prosthetic limb, which optimizes the user's remaining physical capabilities; including range of limb motion and strength. The remaining limb itself, is a very basic limiter of function and performance. Typically, the longer the limb the more functional capacity over a prosthesis will be achieved, assuming the limb's muscular and skeletal framework are not compromised beyond simply hand absence. Loss of muscle tissue, permanent nerve damage and phantom limb sensations-pain can all impact the user's ability to perform with a prosthesis. In general, a user with a limb absence where less than 25% of the humerus remains, will **not** be a good candidate for using upper extremity prosthetic technology successfully. Additionally, the user's cognitive acuity and capability are important in understanding the ranges of function and the specific limits of function, while controlling such prosthetic tools, to optimize their performance and avoid injury.

Finally, all prosthetic terminal devices, have inherent dangers of entanglement or engagement where release can be compromised because of their physical exterior design, unique unforgiving materials, and inanimate lack of "feel". Wearing an upper extremity prosthesis does involve risk! *Training and therapy are always recommended when using a prosthesis, especially when using new or unique technology or changing prosthetic system operative controls.*

Mechanical Lock Gripping Devices:

This class of terminal devices is unique in that they provide the ability to manually lock and unlock around objects, tools and equipment. There is no analog in the human hand for such function. Such function has inherent benefits and risks that the user must anticipate and understand. Certain activities like vehicle control should never be used with devices in this class because instant or quick release of an object or handlebar- steering wheel is not possible. Benefits include very stable control over objects without the liability of releasing or dropping, providing a safety factor for activities where such function is the primary goal.

Storage and Handling

It is recommended that Mechanical Lock Gripping Devices (or prosthetic/orthotic components) are stored in a cool, clean, dry environment away from harsh chemicals (chlorine, acids, acetone, etc.).

Warnings and Precautions

CAUTION: Abnormal or improper environmental conditions will lead to malfunctioning and damage of the prosthesis and is not covered under the warranty of the device. This prosthetic component must not be subjected to dust/debris, liquids other than fresh water, abrasives, vibration, or activities which would damage the biological limb. Do not allow debris or liquids to remain in the prosthesis and its components during use. Rinse the wrist with fresh water and dry immediately after exposure.

CAUTION: These Mechanical Lock Gripping Devices are waterproof to 1 meter; however, if they are submerged, they should be rinsed with fresh water and **dried** immediately to remove salt, chlorine, or debris.

Qualified provider

Attachment, adjustment, alignment, and delivery of this device must be performed by or under the direct supervision of a qualified prosthetist. Unless stated in this manual, any such activities should not be attempted by the user and will potentially void the device warranty.

Specifications and Preparations Before Use

Product Name: BLACK IRON MASTER & TRAINER SERIES

Specifications

	Black Iron Master	Black Iron Trainer	Black Iron Trainer Pivot
Lenght	3.5 in. (9.0 cm.)	2.6 in. (6.6 cm)	3.5 in. (8.9 cm.)
Width	2.0 in (5.0 cm.)	2.0 in. (5.0 cm.)	2.0 in. (5.0 cm.)
Weight	25 oz. (709 gm.)	15.5 oz. (440 gm.)	18 oz. (511 gm.)
Bar Diameters	875 in. (2.2 cm.) - 1.5 in. (3.8 cm.)	875 in. (2.2 cm.) - 1.5 in. (3.8 cm.)	875 in. (2.2 cm.) - 1.5 in. (3.8 cm.)
Load Rating	1650 lbs. (750 kg)	440 lbs. (200 kg)	400 lbs, (200 kg)
Color	Dark Gray/Black/ Silver	Dark Gray/Black/ Silver	Dark Gray/ Black/Silver

Applications

Age:	Adult	Teenage-Adult
	Serious, professional, expert only	All around lifting and training

Limb Description: All Models: Trans-radial & mid-length - long Trans-humeral

PIVOT MODEL: The pivot model provides for a greater range of motion through the exercise cycle. This model is not applicable to wrist disarticulates unless the excess length of the prosthesis is acceptable.

Product Name: BLACK IRON LITE SERIES

Specifications

	BLACK IRON LITE	BLACK IRON LITE PIVOT
Lenght	3.3 in. (8.4 cm.)	4.4 in. (11.0 cm.)
Width	2.0 in. (5.0 cm.)	2.0 in. (5.0 cm.)
Weight	11.7 oz. (332 gm.)	13.0 oz (369 gm)
Bar Diameter	All models: .875 in. (2.2 cm.) - 1.5 in. (3.8 cm.)	
Color	Black, Teal, Violet	

Applications

Age: Both models Teenage-Adult All around lifting and training

Limb Description: All Models: Trans-radial & mid-length - long Trans-humeral

PIVOT MODEL: The pivot model provides for a greater range of motion through the exercise cycle. This model is not applicable to wrist disarticulates unless the excess length of the prosthesis is acceptable.

Product Name: ISHI ARCHERY & FISHING TD

Specifications

	ISHI	ISHI HD
Lenght		3.5 in. (9.0 cm.)
Width		3.1 in. (7.9 cm.)
Weight		1.25 in. (170 gm.)
Opening	Flexible mandibles accommodate a variety of bow and rod diameters	
Color	Stealth Black	

Applications

Age:

ISHI for archery and lighter duty freshwater fishing, 5-10 years old, bows to maximum 40 pounds.

ISHI HD for heavier hunting bows and for heavier duty freshwater and saltwater fishing.

Limb Description: All Models: Trans-radial & mid-length - long Trans-humeral

Product Name: HAMMERHEAD

Specifications

HAMMERHEAD	
Lenght	4.3 in. (11.0 cm.)
Width	2.7 in. (7.0 cm.)
Weight	12.0 oz. (340 gm.)
Opening	Accommodates a wide variety of round and shaped kayak paddle handles
Color	Black

Applications

Age: Teenage through Adult

Limb Description: Trans-radial & mid-length - long Trans-humeral

Product Name: MULTI-D & MULTI-D PIVOT

Specifications

	MULTI-D	MULTI-D PIVOT
Lenght	3.8 in. (9.7 cm.)	4.55 in. (11.6 cm.)
Width	2.0 in. (5.1 cm.)	2.0 in. (5.1 cm.)
Weight	9.6 oz. (272 gm.)	11.4 oz. (323 gm.)
Opening	Adjustable diameters: 1.125 in.(2.86 cm.) - 1.750 in. (4.5 cm) "plus".	
Color	Black	

Applications

Age: Teenage through Adult

Limb Description: All Models: Trans-radial & mid-length - long Trans-humeral

Product Name: FLEXIBLE THING

Specifications

FLEXIBLE THING	
Lenght w/o strap	2.5 in. (6.4 cm.)
Straå Length	12 in. (30.5 cm.)
Width	1.5 in. (3.8 gm.)
Height	1.75 in. (4.4 cm.)
Weight	4.5 oz (128 gm.)
Color	Satin Black

Preparation

Consider the following factors as you determine the most appropriate technology and prosthesis design for your client:

1. Does the user understand the function and safety features of the device?
2. What other terminal devices will be used with this prosthesis?
3. In what types of environments will the device(s) be used?
4. In what positions will the person use the device?
5. How will the user switch terminal devices? Will quick-disconnect be helpful?
6. For what other activities might this device be helpful?
7. Will the person benefit from additional training in the use of the device?

Installation

These Mechanical Lock Gripping Devices may be installed in any Fillauer TRS Omega Wrist, or any Fillauer LLC wrist unit with a ½-20 thread. Follow the instructions provided with the wrist unit for best results.

Cabling

These Mechanical Lock Gripping Devices are designed to be used without any prosthetic cable systems.

Compatibility

Fillauer TRS, Mechanical Lock Gripping Devices have been evaluated with and are recommended for use with Fillauer TRS Omega Wrists, and Fillauer LLC wrists that have a ½-20 internal thread. They may be used with any equivalent ½-20 threaded wrist units; however, damage caused by other manufacturers wrist units is not covered under warranty of this device.

Care and Maintenance

General:

Upper Extremity Terminal Devices-End Effectors:

Prosthetic devices should be used and treated as quality tools, regularly inspected for damage, deterioration and or excessive wear. The devices are waterproof and or highly water resistant and constructed of materials that will not rust or corrode or are resistant to such exposure. They can be immersed in water and washed. If truly mechanical they should be dried, as thoroughly as possible, and light lubricants can be applied to moving mechanical linkages as needed. Soap and water cleansing or rinsing are recommended, especially after operation in environments with a high quantity of dirt, dust, grit, oils, or solvents. If used in ocean water, cleaning as soon as possible in fresh water is highly recommended with all devices. Long term exposure in intense direct sunlight or heat outside of normal use will damage prosthetic devices that incorporate synthetic rubber in their construction or coverings.

Mechanical Lock Gripping Devices:

These devices employ simple reliable locking type mechanisms. With extended or heavy use certain parts can wear or deteriorate, so inspect the devices regularly before each use. Clean as needed. These devices require little if any maintenance.

Disposal/Waste Handling

The product must be disposed of in accordance with applicable local laws and regulations. If the product has been exposed to bacteria or other infectious agents, it must be disposed of in accordance with applicable laws and regulations for the handling of contaminated material.

All metal components may be removed and recycled at the appropriate recycling facility.

User Information

The providing health care professional must review the following information directly with the user:

Consider the following factors as you prepare for this device:

1. Are you able to don and to doff the device independently?
2. Are you able to attach/detach the cable from the device?
3. Are you able to switch terminal devices independently?
4. Are you able to access services from a skilled occupational therapist for additional training if needed?

Warnings and Precautions for the user

CAUTION: Abnormal or improper environmental conditions will lead to malfunctioning and damage of the prosthesis and is not covered under the warranty of the device. This prosthetic component must not be subjected to dust/debris, liquids other than fresh water, abrasives, vibration, or activities which would damage the biological limb. Do not allow debris or liquids to remain in the prosthesis and its components during use. Rinse the wrist with fresh water and dry immediately after exposure.

CAUTION: These Mechanical Lock Gripping Devices are waterproof to 1 meter; however, if they are submerged, they should be rinsed with fresh water and dried immediately to remove salt, chlorine, or debris.

Serious Incidents

In the unlikely event of a failure resulting in a fall and/or injury, seek immediate medical help and contact your prosthetist at the earliest possible convenience

Warranty

This product has a 12-month warranty against manufacturer defects

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© 2021 Motion Control, Inc. 1910013 Rev H 06-29-2021
Motion Control, Inc.