

Economy Friction Wrists

Product Manual

The Fillaauer logo is written in a white, cursive script font on a blue rectangular background. To the left of the logo, there is a faint, light blue line drawing of a wrist device with concentric rings and a textured surface.

Intended Use

Economy Friction Wrists have variable friction through the use of compressible discs. Increased compression increases the resistance to rotation. These wrists are not recommended for terminal devices with large round wrist bases.* The large diameter base can prevent the terminal device's threaded adapter from penetrating deep enough through the rubber friction element to engage the threaded receiver effectively. The result of such an installation is that the threads of the terminal device do not fully engage into the wrist, thereby compromising the stability of the TD and in most cases causing premature thread wear or breakage.

***JAWS, Grip V/C TDs, ADEPT TDs, Lite-Touch TDs, and many activity specific TDs such as Hammerhead, BIT, Swinger, Multi-D, etc.**

Performance Characteristics

- Available with aluminum, stainless steel, or titanium body
- Stainless steel-threaded bushing fits into the body
- Titanium version is lightweight and corrosion resistant
- ½-20 thread

Economy Wrist

Material	Diameter	Build Height	Weight	Lamination Ring	Product Number
Al	2.0 in. (51 mm)	1.2 in. (30 mm)	2.3 oz. (64 g)	No	52159
SS	2.0 in. (51 mm)	1.2 in. (30 mm)	4.7 oz. (134 g)	No	52162
Ti	2.0 in. (51 mm)	1.2 in. (30 mm)	3.3 oz. (93 g)	No	59491
Al	2.0 in. (51 mm)	1.26 in. (32 mm)	2.4 oz. (68 g)	Yes	58690

Oval Economy Wrist

Material	Length × Width	Build Height	Weight	Lamination Ring	Product Number
Al	2 ¼ × 1 ½ in. (56 × 41 mm)	1.38 in. (35 mm)	2.5 oz. (72 g)	No	51622
Al	2 ¼ × 1 ½ in. (56 × 41 mm)	1.24 in. (31 mm)	2.7 oz. (76 g)	Yes	58700

Economy Wrist with Straps

Material	Diameter	Build Height	Weight	Product Number
SS	2.0 in. (51 mm)	1.2 in. (30 mm)	6.0 oz. (170 g)	52163
Ti	2.0 in. (51 mm)	1.2 in. (30 mm)	4.9 oz. (140 g)	59491

The device is intended for single user/patient use only.

Storage and Handling

It is recommended that prosthetic wrists be stored in a cool, clean, dry environment away from harsh chemicals (chlorine, acids, acetone, etc.).

Warnings and Precautions



NOTICE: An upper-limb prosthetic device user's ability to drive should be determined on a case-by-case basis by a specialist. Contact your local governing authorities regarding any driving restrictions or limitations.



WARNING: Body-powered devices should not rely on cable tension for grasp control if the user has been cleared to drive with the prosthesis. Failure to maintain tension while controlling the steering wheel could cause serious injury or death.



CAUTION: Do not tighten the friction screw without first inserting a terminal device with a 1/2-20 thread.



CAUTION: Abnormal or improper environmental conditions will lead to malfunctioning and damage of the prosthesis and are not covered under the warranty of the device. This prosthetic component must not be subjected to dust/debris, liquids other than fresh water, abrasives, vibration, activities which would damage the biological limb, or prolonged extreme temperatures (< -5 °C or > 50 °C). Do not allow debris or liquids to remain in the prosthesis and its components during use. Rinse the wrist with fresh water and dry immediately after exposure.



CAUTION: The wrist unit is waterproof to 1 meter. However, if the wrist is submerged, it should be rinsed with fresh water and dried immediately to remove salt, chlorine, or debris.

Qualified Provider

Attachment, adjustment, alignment, and delivery of this device must be performed by or under the direct supervision of a qualified prosthetist. Unless stated in this manual, any such activities should not be attempted by the user and will potentially void the device warranty.

Specifications and Preparations Before Use (Risk Management for Installation and Calibration)

Alignment

Prosthetic wrists should be aligned to provide the best possible work envelope for the patient's specific goals. Standard alignment begins at 5° of flexion and 5° of radial deviation but should be tailored to the individual patient.

Fabrication Without Lamination Ring

1. Check to make sure that the wrist can be placed on the distal end of the socket in such a way that the internal components do not contact the inner socket before the wrist body.
2. The wrist may be laminated by removing the cover and internal components from the wrist body.
3. Place the wrist body on the distal end of the forearm mold (beeswax, foam, plaster, or similar).
4. A PVA bag should be used to separate the wrist body from the forearm mold if foam or plaster are used.
5. Wax the interior surface of the wrist body and the distal face of the body.
6. Pack the wrist body with silicone fitting gel or similar to prevent it from filling with laminate.
7. Mask the body distal to the tie-in groove on the lamination ring to preserve access to any moving components and to keep all laminate clear of the distal end.
8. Laminate with the appropriate materials for durability and finish as desired by the patient, being sure to tie each structural layer into the tie-in groove in the lamination ring.
9. Carbon fiber tape is a good choice for reinforcing the connection to the lamination ring and is commonly used as distal to proximal strips tied in with circumferential wraps.

Fabrication With Lamination Ring

1. Place the wrist unit on the distal end of the forearm mold (beeswax, foam, plaster, or similar). Ensure that the internal components do not make contact with the inner socket.
2. The wrist may be laminated by removing the face assembly and internal components from the lamination ring once its alignment and clearance are confirmed.
3. A PVA bag should be used to separate the lamination ring from the forearm mold if foam or plaster are used.
4. Wax the interior surface of the lamination ring.
5. Pack the lamination ring with silicone fitting gel or similar to prevent it from filling with laminate.
6. Mask the lamination ring distal to the tie-in groove to keep all laminate clear of the distal end and away from any moving components.
7. Laminate with the appropriate materials for durability and finish as desired by the patient, being sure to tie each structural layer into the tie-in groove in the lamination ring.
8. Carbon fiber tape is a good choice for reinforcing the connection to the lamination ring and is commonly used as distal to proximal strips tied in with circumferential wraps.
9. When the lamination is complete clean the lamination ring and remove any packing material. Ensure that no lamination is left distal to the tie-in groove.
10. Slide the nutplate through laminating ring and pull it into place. Keep the engaging bumps proximal.
11. When reattaching the face assembly and/or internal components, be sure to use a medium-strength threadlocker on the attachment screws before delivery.
12. Install the face assembly, engaging the nutplate with the #4-40 screws. Tighten the screws.
13. If the orientation of the button/lever needs to be changed, loosen the #4-40 screws approximately 1/16 inch (2 mm) and

gently tap the face assembly. To loosen the nutplate, rotate to the proper orientation and re-tighten.

Adjusting the Lamination Ring


1. If the wrist unit has a separate lamination ring, its alignment can be adjusted during the installation of the internal components.
2. To install, slide the nutplate through laminating ring and pull it into place. Keep the engaging bumps proximal.
3. When reattaching the face assembly and/or internal components, be sure to use a medium-strength threadlocker on the attachment screws before delivery.

4. Install the face assembly, engaging the nutplate with the #4-40 screws. Tighten the screws.
5. If the orientation of the button/lever needs to be changed, loosen the #4-40 screws approximately $\frac{1}{16}$ inch (2 mm) and gently tap the face assembly. To loosen the nutplate, rotate to the proper orientation and re-tighten.

Wrist Assemblies


52159 WE-500 Aluminum
52162 WE-500S Stainless Steel
59491 Titanium

- A.** 52164 Body, WE-500
 52831 Body, WE-500S
 59490 Body, Titanium Wrist
- B.** 52171 Rubber
- C.** 52172 Cover, WE-500
 59464 Cover, WE-500S, and Titanium Wrist
- D.** 52173 Brass Washers




61622 OW-100 Adult Size Aluminum

- A.** 51624 Body, 51622
- B.** 52171 Rubber
- C.** 52172 Cover
- D.** 52173 Brass Washer




58690 WE-500LR, Aluminum Wrist with Laminating Ring

- A.** 59522 Laminating Ring
- B.** 58693 Body, 58690
- C.** 58692 Nut Plate
- D.** 58283 Screw, 4 Each
- E.** 52171 Rubber
- F.** 52172 Cover
- G.** 52173 Washer, 2 Each



58700 OW-100LR with Laminating Ring

- A.** 59526 Laminating Ring
- B.** 58698 Body Assembly
- C.** 52171 Bushing
- D.** 52172 Cover
- E.** 52173 Washer, 2 Each
- F.** 70242 Screw, #8-32, 2 Each



Disposal / Waste Handling

The product must be disposed of in accordance with applicable local laws and regulations. If the product has been exposed to bacteria or other infectious agents, it must be disposed of in accordance with applicable laws and regulations for the handling of contaminated material.

All metal components may be removed and recycled at the appropriate recycling facility.

Warranty

This product has a 12-month warranty against manufacturer defects.

User Instructions

The providing health care professional must review the following information directly with the user.

Warnings and Precautions for the User



NOTICE: The user should monitor their prosthesis daily and contact their health care professional if they experience changes in device performance or if it begins to make noise



CAUTION: All maintenance should be performed by the qualified health care professional.



NOTICE: An upper-limb prosthetic device user's ability to drive should be determined on a case-by-case basis by a specialist. Contact your local governing authorities regarding any driving restrictions or limitations.



CAUTION: Body-powered devices should not rely on cable tension for grasp control if the user has been cleared to drive with the prosthesis. Failure to maintain tension while controlling the steering wheel could cause serious injury or death.



CAUTION: Do not tighten the friction set screw without first inserting a terminal device with a 1/2-20 thread.



CAUTION: Abnormal or improper environmental conditions will lead to malfunctioning and damage of the prosthesis and are not covered under the warranty of the device. This prosthetic component must not be subjected to dust/debris, liquids other than fresh water, abrasives, vibration, activities which would damage the biological limb, or extreme temperatures (< -5 °C or > 50 °C). Do not allow debris or liquids to remain in the prosthesis and its components during use. Rinse the device with fresh water and dry immediately after exposure.



CAUTION: The wrist unit is waterproof to 1 meter. However, if the wrist is submerged, it should be rinsed with fresh water and dried immediately to remove salt, chlorine, or debris.

Serious Incidents

In the unlikely event a serious incident occurs in relation to the use of the device, users should seek immediate medical help and contact their prosthetist, local competent authority, and Fillauer at the earliest possible convenience. Clinicians should at any time contact their local Fillauer representative and local competent authority immediately in the event of any device failure.



To see other language options, visit fillauer.com.

Fillauer®

www.fillauer.com



Fillauer LLC

2710 Amnicola Highway
Chattanooga, TN 37406
423.624.0946



Fillauer Europe

Kung Hans väg 2
192 68 Sollentuna, Sweden
+46 (0)8 505 332 00



© 2022 Fillauer LLC
M097/06-08-22/Rev.1