
MC ProWrist Rotator

Protetikermanual

Fillauer[®]
Motion Control

Særlige forholdsregler

Risikostyring

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
Protetikermanual

Introduktion

The Motion Control ProWrist (Figure 1) provides powered wrist rotation with an on-board microprocessor. This controller provides Proportional Control and switching of function between the terminal device and wrist rotation. The ProWrist is to be used in cases of transradial and higher-level amputation. The ProWrist may also be used with body-powered elbows and some other manufacturers' externally-powered elbows for transhumeral and higher-level amputations. When using a powered elbow, the MC Standard version of the wrist rotator should be used.

Equipped with an industry-standard quick disconnect wrist unit and coaxial plug, the ProWrist is compatible with almost all available powered terminal devices. The optional Wrist Disable Switch (p/n 3010846, Figure 3) allows disabling of Hand/Wrist switching for times when the patient desires only the terminal device to be operational.

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Systemkrav

- iPad, iPhone, or iPod Touch running iOS 11 or later.

Download the
MAUI software on
the Apple Store.

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Bestillingsoplysninger

Description	Tan	Brown	Jet Black
MC Motor Direct Wrist Rotator, 4-Band Sequential			5010045
MC Motor Direct Wrist Rotator, 6-Band Sequential			5010159
MC ProWrist Rotator, 4-Band Sequential — 2-channel			5010056
MC ProWrist Rotator, 6-Band Sequential — 2-channel			5010157
MC ProWrist Rotator, 6-Band Simultaneous — 4-channel			5010158
ProWrist Disable Switch			3010846
Lamination Collar for MC Wrist	1100292	1100296	1100288
Lamination Collar for MC Wrist for Hosmer Prefab (Large), Boston or AFB	1100295	1100299	1100291
Forearm Endcap			3010758

Use a Motor Direct Wrist Rotator for these elbows:

Utah Arm (U3/U3+), E2 Elbow with pattern recognition, Ottobock Dynamic Arms, Steeper Espire, other manufacturer's powered elbows. Devices using Coapt and Glide Systems in most configurations.

Use a ProWrist Rotator for these elbows:

Motion Arm Hybrid, E2 Elbow standard inputs. Ottobock Ergo, Steeper Classic, etc.

4-Band Terminal Devices:

TD's running open/close dual site systems, ilimb Quantum, Adams Hand, Vincent, bebionic, MC Motor Direct and ProPlus TD's, other manufacturer TD's.

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Statisk belastning: 22,7 kg. (50 lbs)

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WHICH WRIST ROTATOR DO I NEED?

Wrist Configuration	Coaxial Plug	Has Processor	Usage	Required Items	Compatible Elbow Systems and Devices
Motor Direct Part: 5010045	4-band Sequential	NO	Pattern Recognition, powered elbows, devices with other upstream configurations. 4-band used with most dual site systems, any device that operates open/close.	50mm Lamination collar: Smooth (device) or Slotted (socket). Lamination Dummy.	Coapt, Glide, Utah Arm, E2 Elbow pattern system, Ottobock Dynamic Arm or Dynamic Arm Plus, Steeper Espire Pro
Motor Direct 6-Band Part: 5010159	6-Band Sequential	NO	Pattern Recognition, feedback devices, powered elbows, devices with other upstream configurations.	50mm Lamination collar: Smooth (device) or Slotted (socket). Lamination Dummy.	Coapt, Glide, Utah Arm, E2 Elbow pattern system, Ottobock Dynamic Arm or Dynamic Arm Plus, Steeper Espire Pro.
ProWrist Part: 5010056	4-band Sequential, 2 channel	YES	Dual Site TransRadial System. 4-band used with most dual site systems.	50mm Lamination collar: Smooth (device) or Slotted (socket). Lamination Dummy. Optional: ProWrist Disable Switch.	Motion Arm, E2 Elbow standard input, hybrid elbows, Ottobock Ergoarm, Steeper Espire Hybrid, Classic and Classic Plus.
ProWrist 6-Band Seq. Part: 5010157	6-Band Sequential, 2 channel	YES	Dual Site TransRadial System. 6-band needed for vibration feedback.	50mm Lamination collar: Smooth (device) or Slotted (socket). Lamination Dummy. Optional: ProWrist Disable Switch.	Motion Arm, E2 Elbow standard input, hybrid elbows, Ottobock Ergoarm, Steeper Espire Hybrid, Classic and Classic Plus.
ProWrist 6-Band Sim. Part: 5010158	6-Band Simultaneous, 4 channel	YES	Wearer uses 4 electrodes (two sets) or implantable electrodes (not common use). Must use multiple inputs. 6-band needed for vibration feedback.	50mm Lamination collar: Smooth (device) or Slotted (socket). Lamination Dummy.	Motion Arm, hybrid elbows, Ottobock Ergoarm Hybrid Plus and Electronic Plus (50mm), Steeper Espire Hybrid, Classic and Classic Plus.

Customer Support

Americas, Oceania, Japan

Fillauer Motion Control

115 N. Wright Brothers Drive
Salt Lake City, UT 84116
+1 801-326-3434
motioninfo@fillauer.com

Europe, Africa, Asia

Fillauer Europe

Kung Hans väg 2
192 68 Sollentuna, Sweden
+46 (0)8 505 332 00
support@fillauer.com

www.fillauer.com

Fillauer

2710 Amnicola Highway
Chattanooga, TN 37406
+1 423.624.0946

Fillauer Europe

Kung Hans väg 2
192 68 Sollentuna,
Sweden
+46 (0)8 505 332 00