

Wrists

Economy Wrists

Economy Wrist	Product Number	Diameter
	52159 Al 64g	2 in. 5.1 cm
	52162 SS 134g	2 in. 5.1 cm
OW Oval Wrist, Economy	Product Number	Length × Width
	58700 Al 76g LR	2 ¼ × 1 ⅝ in. 5.7 × 4.1 cm
WA Disarticulation Friction Wrist	Product Number	Diameter
	52119 SS 57g	1 ¾ in. 3.5 cm

Constant Friction Wrists

WE Infant Friction Wrist	Product Number	Diameter
	52142 Al 28g NR	1 ¼ in. 3.2 cm
WE Friction Wrist	Product Number	Diameter
	52147 Al 42g	1 ½ in. 3.8 cm
	52156 Al 57g	1 ¾ in. 4.4 cm
	52160 Al 74g	2 in. 5.1 cm
OW Oval Friction Wrist	Product Number	Length × Width
	51618 Al 57g	2 × 1 ½ in. 5.1 × 3.8 cm
	51623 Al 71g	2 ¼ × 1 ⅝ in. 5.7 × 4.1 cm
OW Friction Wrist Insert	Product Number	Dimensions
	52151 Al 30g	½-20 Insert

SS Stainless Steel	ABS ABS
Al Aluminum	g Product Weight, includes Insert
Ti Titanium	LR Lamination Ring Model
D Delrin	NR Non-removable Insert

Shapeable Passive Hand Wrist	Product Number	Length × Width
	62021 ABS 28g	2.2 × 1.5 in. 5.6 × 3.8 cm
CAPP Wrist	Product Number	Diameter
	55108 D 16g	1 ½ in. 3.8 cm
	55109 D 22g	1 ¾ in. 4.4 cm
	55110 D 30g	2 in. 5.1 cm

Flexion Wrists

APRL Child Flexion Wrist	Product Number	Diameter
	52176 Al 42g	1 ¼ in. 3.2 cm
FW Flexion Friction Wrist	Product Number	Diameter
	51084 Al 71g	1 ¾ in. 4.4 cm
	51100 Al 92g	2 in. 5.1 cm
Sierra Flexion Wrist	Product Number	Diameter
	52371 Al 113g	2 in. 5.1 cm

Quick Disconnect Locking Wrists

Quick Disconnect Wrist	Product Number	Diameter
	62342 Al 89g	1 ¾ in. 4.4 cm
	61862 Al 105g	2 in. 5.1 cm
Quick Disconnect Oval Wrist	Product Number	Length × Width
	62343 Al 92g	2 × 1 ½ in. 5.1 × 3.8 cm
Quick Disconnect Insert	Product Number	Dimensions
	61860 SS 23g	½-20 Insert
FM Quick Change Wrist	Product Number	Diameter
	51043 Al 92g	2 in. (5.1 cm)
	51044 SS 170g	2 in. (5.1 cm)
	58297 Al 112g LR	2 in. (5.1 cm)
	62436 Ti 166g LR	2 in. (5.1 cm)
FM Quick Change Wrist Insert	Product Number	Dimensions
	51052 SS 29g	½-20 Insert

OMEGA Wrists

OMEGA Wrist	Product Number	Diameter
	● OMEGA R ● OMEGA RB Al-SS 103g	2 in. 5.2 cm
OMEGA Oval Wrist	Product Number	Length × Width
	● OMEGA M ● OMEGA MB Al-SS 89g	2 × 1 ½ in. 5.1 × 3.8 cm
	● OMEGA L ● OMEGA LB Al-SS 111g	2.4 × 1.7 in. 6.1 × 4.3 cm
OMEGA Insert	Product Number	Dimensions
	OM003 SS 23g	½-20 Insert