



Xtra Torsion. Xtra Shock.

AllPro XTS Foot

**Foot
Weight**
25.6 / 27.8 oz.
728 / 789 g

**Maximum
Patient Weight**
276 lbs.
125 kg

**Build
Height**
6.4 / 8.4 in.
16.2 / 21.3 cm

**Activity
Level**
K3 / K4

Fillauer®



More Torsion in a Lower Profile

It's the AllPro performance you have come to expect with an even smoother rollover and forgiveness on rough terrain or around town.

The XTS unit provides adjustable torsion and vertical shock reduction without sacrificing build height. This field-serviceable unit is built with an active lifestyle in mind.



Waterproof

Have the confidence to go where you want. The AllPro XTS is IP67 waterproof from tip to toe.



Industry-Leading Flexibility


The split-toe, carbon composite design provides multiaxial movement and dynamic response.



AllPro XTS

Specifications

	6.4-inch model	8.4-inch model
Patient Weight	275 lbs. (125 kg)	275 lbs. (125 kg)
Foot Build Height	6.4 in. (16.3 cm)	8.4 in. (21.3 cm)
Foot Weight (26 cm, Category C)	25.6 oz. (728 g)	27.8 oz. (789 g)
Functional Level	K3 – K4	
Primary Materials	Carbon Composite, Aluminum, Polyurethane, and Stainless Steel	
Durable	Meets ISO-22675 Standard	
Waterproof	The foot unit is waterproof to 1 meter.	
Warranty	36 months	
Patient Trial	60 days	



Ordering Information

Foot Shell sold separately.

- 229-106-XX-XX** AllPro XTS foot, 6.4 in. pylon, foot sizes 22 – 30 cm
- 229-108-XX-XX** AllPro XTS foot, 8.4 in. pylon, foot sizes 22 – 30 cm

Patient Weight		22 cm	23 cm	24 cm	25 cm	26 cm	27 cm	28 cm	29 cm	30 cm
100 – 119 lbs.	45 – 54 kg	22-A2	23-A2	24-A2	25-B4	26-B4				
120 – 139 lbs.	54 – 63 kg	22-A3	23-A3	24-A3	25-B4	26-B4	27-C6	28-C6		
140 – 159 lbs.	64 – 72 kg	22-B4	23-B4	24-B4	25-B4	26-B4	27-C6	28-C6	29-D7	30-D7
160 – 179 lbs.	73 – 81 kg	22-B4	23-B4	24-B4	25-C5	26-C5	27-C6	28-C6	29-D7	30-D7
180 – 199 lbs.	82 – 90 kg	22-C5	23-C5	24-C5	25-C5	26-C5	27-D7	28-D7	29-D7	30-D7
200 – 224 lbs.	91 – 102 kg	22-C5	23-C5	24-C5	25-D6	26-D6	27-D7	28-D7	29-E8	30-E8
225 – 250 lbs.	102 – 113 kg				25-D6	26-D6	27-E8	28-E8	29-E8	30-E8
250 – 276 lbs.	113 – 125 kg				25-E7	26-E7	27-E8	28-E8	29-F9	30-F9

AllPro XTS Cosmetic Cover Options

ProCover



Split Design

ProCover is easy to install and remove for cleaning and maintenance.



Innovative Materials

The durable polyurethane foam maintains superior energy return over time — ideal for high activity patients.



Design-to-Fit Shape

The unique shape enhances rollover performance and provides a secure fit for a more responsive feel.

Specifications

Product Weight (26 cm)	Functional Level	Material	Waterproof	Warranty
9.1 oz. (259 g)	K3 – K4	Polyurethane foam	Waterproof to 1 meter	9 months

Ordering Information

ProCover purchase includes two hardware kits.

	22 cm	23 cm	24 cm	25 cm	26 cm	27 cm	28 cm	29 cm	30 cm
ProCover	182-10-1022	182-10-1023	182-10-1024	182-10-1025	182-10-1026	182-10-1027	182-10-1028	182-10-1029	182-10-1030

Replacement Parts

	Quantity	Size	Part Number
Hardware Replacement Kit	4 Connector Pins, 2 Velcro® Straps	22 – 25 cm	182-10-2300
Hardware Replacement Kit	4 Connector Pins, 2 Velcro® Straps	26 – 30 cm	182-10-2350

Patent US10874529B2.

Foot Shell

Made with BIOVYN™, a bio-attributed PVC that is based on 100% renewable feedstock, our product not only contributes to significantly reduced carbon emissions but also offers enhanced durability. To complement the foot shell, a sleek, black prosthetic sock, sized appropriately, is now included.



Specifications

Product Weight (26 cm)	Material	Warranty
8.9 oz. (251 g)	BIOVYN™ PVC	9 months

Ordering Information

Replace **XX** with foot length (22 – 30 cm). Foot shell includes black prosthetic sock.

Color	Left	Right
03	45- XX -13-03-3	46- XX -13-03-3
09	45- XX -13-09-3	46- XX -13-09-3
13	45- XX -13-13-3	46- XX -13-13-3
99	45- XX -13-99-3	46- XX -13-99-3

Prosthetic Socks

Size	Part Number
22 – 24 cm	636202224
25 – 27 cm	636202527
28 – 30 cm	636202830

Fillauer®

Fillauer LLC

A Fillauer Company

2710 Amnicola Highway
Chattanooga, TN 37406
423.624.0946

© 2023 Fillauer LLC
AD392A/11-23-22/08-25-23