

**CAUTION - Be sure to use the proper size pyramids and receivers together.**

**Pyramid sizes are defined as:**

- **SMALL** Mates with 20mm spherical radius (commonly used throughout many pediatric systems including on top of the knee).
- **MEDIUM** Mates with 22mm spherical radius (standard pyramid for MightyMite system, also used at the ankle in some pediatric systems).
- **LARGE** Mates with 30mm spherical radius (commonly used throughout adult systems).

**MightyMite pyramid and receiver components are compatible with:**

FC100110	SACH Foot Adapter w/Med Pyramid	Medium Pyramid
FC100114	Four Hole Adapter w/Med Pyramid	Medium Pyramid
FC100118	Integrated Pylon w/Aluminum Receiver	Medium Pyramid
FC100119	Integrated Steel Pylon w/Titanium Receiver	Medium Pyramid
FC100120	Integrated Steel Pylon w/Aluminum Receiver	Medium Pyramid
FC100121	Integrated Steel Pylon w/Aluminum Receiver	Large Pyramid
FC100122	Aluminum Tube Clamp f/Medium Pyramid	Medium Pyramid
FC100123	Titanium Tube Clamp f/Medium Pyramid	Medium Pyramid
FC100124	Aluminum Tube Clamp f/Large Pyramid	Large Pyramid
FC100134	Rotatable Medium Pyramid Receiver	Medium Pyramid
FC100144	Rotatable Small Pyramid Receiver	Small Pyramid
FC100152	BK Socket Adapter	Medium Pyramid
FC100530	M30 Thread to Pyramid Adapter	Medium Pyramid



Fillauer LLC • 2710 Amnicola Highway • Chattanooga, TN 37406 • Tel: 800.251.6398 • [www.fillauer.com](http://www.fillauer.com)  
Hosmer • 561 Division Street • Campbell, CA 95008 • Tel: 800.827.0070 • [www.hosmer.com](http://www.hosmer.com)  
Centri® • Kung Hans Väg 2 • 192 68 Sollentuna, Sweden • Tel: +46 8 505 332 00 • [www.centri.se](http://www.centri.se)

PM195/04-09/05-10

