



PARTS CATALOG

UTAH ARM 3

UTAH ARM 3 PACKAGE Tan: 5010035 Brown: 5010036 Jet Black: 5010038
UTAH ARM 3+ PACKAGE Tan: 5010039 Brown: 5010040 Jet Black: 5010041




EXTENSION CABLE (EMG)
Short 10" (25cm): 3010094
Long 6' (1.5m): 3010073




RIBBON CABLE U3 CONNECTOR
Single Channel: 300753
Dual Channel: 3010752
9-Pin Single: 3010731


ALTERNATE INPUTS:




LINEAR POTENTIOMETER
3010546



TOUCH PAD
Tan: 3010669
Brown: 3010670
Jet Black: 3010721




1-CHANNEL TOUCH PAD CABLE
3010548




2-CHANNEL TOUCH PAD CABLE
3010549


WIRE HARNESS:




1-WIRE HARNESS
3010544
• Hand & Elbow



3-WIRE HARNESS
3010540
• Hand
• Elbow
• Wrist Switch



5-WIRE HARNESS
3010545
• Hand & Elbow
• Wrist Switch
• Remote Lock/Unlock
• Remote Power Switch




BATTERY CHARGER
3010732

POWER SUPPLY
1390040


POWER CORDS
US: 1390043, EU: 1390043

SINGLE BATTERY (LI-ION)
Tan: 3010745, Brown: 3010746,
Jet Black: 3010747


SUPPLEMENTAL BATTERY PACK/CHARGER
(Charger, Power Supply, Cord, 2 Batteries)
Tan: 3010856, Brown: 3010857,
Jet Black: 3010858



LAMINATION COLLAR
Aluminum: 1070004,
Temp. Plastic: 2000139



FRICTION BAND (Single)
1100015



HUMERAL CLAMP BAND ASSEMBLY
(Friction Band and Clamp Band included)
Tan: 3010145, Brown: 3010146, Jet Black: 3010700

HUMERAL ROTATOR ASSEMBLY
(Lamination Collar, Friction and Clamp Bands included)
Tan: 3010035, Brown: 3010147, Jet Black: 3010701


HUMERAL PLASTIC LAMINATION COLLAR ALIGNMENT FIXTURE
(Set) 3010714




BODY POWERED ELBOW ADAPTER PLATE
1070076 - Adapts lamination collar and
clamp band to E400 Hosmer Elbow

MYO/ONE ELECTRODES:

Standard 10" Cable
Long 18" Cable



U3




3-Pin Available

SNAP U3 STANDARD 5010155 (3mm post), 5010156 (6mm post)
CONTROL SEAL® U3 CABLE STANDARD 5010153
CONTROL SEAL® U3 CABLE LONG 5010154
Single electrode sets and 3-Pin versions of the Myo/One
are also available, see website.


"A" WRIST CONTROL (12-volt proportional control)



MC PULL SWITCH Jet Black: 3010259



HIGH PROFILE WRIST PUSH SWITCH
T: 3010250, B: 3010251, JB: 3010715




LOW PROFILE WRIST PUSH SWITCH
T: 3010675, B: 3010676, JB: 3010720

**ADAPTER CABLE FOR "A" OTTOBOCK
WRIST SWITCH** 3010105
EXTENSION CABLE 9"(23cm) 3010556

LOCK/UNLOCK SWITCHES



MC PULL SWITCH Jet Black: 3010261



HIGH PROFILE PUSH SWITCH
T: 3010252, B: 3010253, JB: 3010717




LOW PROFILE PUSH SWITCH
T: 3010673, B: 3010674, JB: 3010722

ADAPTER CABLE FOR OTTOBOCK SWITCH
3010107

REMOTE POWER SWITCHES



REMOTE PULL SWITCH JB: 3010263




HIGH PROFILE REMOTE PUSH SWITCH
T: 3010254, B: 3010255 JB: 3010716

LOW PROFILE REMOTE PUSH SWITCH
T: 3010671, B: 3010672, JB: 3010723

TRIAD PREAMPS:



SINGLE PREAMP
3010754



ELECTRODE SET (6)
Short: 3010636
Med: 3010835
Long: 3010834



PREAMP DUMMY
1100278



PREAMP GROMMET
1100272



DUMMY PREAMP GROMMET
1100279

SETS:


DUAL CHANNEL UTAH ARM CONNECTOR
3010782

SINGLE CHANNEL UTAH ARM CONNECTOR
3010783


TRIAD SET-NO CABLE
3010786
Single electrodes available see website.

***COAPT® or IBT PATTERN RECOGNITION**
If the system is going to use COAPT or IBT,
please note that during the order.


"B" WRIST CONTROL
(6-volt on/off control) Use with Wrist Control
Adapter Cable for Version B (3010106)




CABLE PULL SWITCH
Ottobock 1701023



HARNESS PULL SWITCH
Ottobock 1701022



ROCKER SWITCH
Ottobock 1701021



**ADAPTER CABLE FOR
"B" OTTOBOCK WRIST
SWITCH** 3010106

EXTENSION CABLES

HALF-MOON EXTENSION CABLES


	2-Pin	3-Pin
12" (30.5cm):	3010951	3010989
18" (46cm):	3010842	3010990
24" (60cm):	3010952	3010991

4-PIN EXTENSION CABLES

12" (30.5cm):	3010325
24" (61cm):	3010326
36" (91.5cm):	3010304

**9-PIN PREAMP
EXTENSION**
3010336

WRIST ROTATORS



MOTOR DIRECT
4-Band: 5010045
6-Band: 3010869

PROWRIST
4-Band: 5010056
6-Band Sequential: 3010832
6-Band Simultaneous: 3010833

PROWRIST DISABLE SWITCH
3010846

WRIST SWITCH EXTENSION CABLE
3010556

WRIST LAMINATION DUMMY
1070495

U3 UPGRADES

**MULTI-ARTICULATING HAND SUPPLEMENT
BATTERY UPGRADE:** 3010677

COAPT® PATTERN RECOGNITION UPGRADE
3010891

WRIST AND TD OPTIONS AVAILABLE
ON WWW.FILLAUER.COM