


Tools, Equipment & Materials





In 1914, George Fillauer Sr. founded the Red Star Pharmacy in Chattanooga, Tennessee. In 1934, his son and fellow CPO, Carlton Fillauer expanded his father's business to establish Fillauer Orthopedic as an innovator in the orthotics and prosthetics industry. Today, the Fillauer Companies are headed by Karl Fillauer as Chairman and his son Michael as the CEO, since their inception, the Fillauer Companies have always been headed by CPO's, and have grown to be a globally renowned leader in the development, manufacturing, distribution and application of orthotic and prosthetic products. The expansion has mainly taken place by acquisitions of smaller manufacturing companies with a unique profile. Over the years we have become nine companies under the Fillauer umbrella and during the 100 year celebration in 2014 a decision was taken on how the company would be steered into the next century. Therefore a consolidation took place in 2015 where all the companies became Fillauer.

Fillauer Europe was founded as Centri in 1948. Located just outside Stockholm, Sweden, Centri joined the Fillauer Companies 2002. While Centri maintained the close, cooperative relationships it traditionally had with prosthetists in Scandinavia, as Fillauer Europe it has now become the European headquarters of Fillauer. With its extensive portfolio of products, Fillauer Europe is committed to working with prosthetists both in Europe, but also in other regions of the world in order to create only the best solutions for patients with all manner of prosthetic and orthotic needs. In addition, Fillauer Europe plays a lead Research & Development role within Fillauer. Fillauer Europe's products include upper and lower extremity prosthetics and the Dynamic Walk which has catapulted us to the forefront of the drop foot orthotics market. Our R&D department also works in partnership with many of Sweden's leading hospitals to develop innovative and high-tech machinery.



Index

Catalog symbol legend	6
-----------------------	---

Trautman Carvers & Sanders

Trautman Double Drum Sander	9
Trautman Variable Speed Floor Carver	10
Trautman Single Speed Floor Carver	11
Trautman Tech Station	12
Dust Collector for Floor Carver	13
Dust Collection Hose Kit for 2100L	13
Air Attachment for Floor Carver	14
Accessory Holder for Floor Carver	14
Lift Assist Kit for Floor Carver	15
Trautman Variable Speed Bench Carver	16
Trautman Single Speed Bench Carver	16
Bench Carver Floor Stand	17
Accessory Holder for Bench Carver	17
Dust Collector for Bench Carver	17
Drill Attachments 1/2"-13 Thread	18
Soft Foam Shaper	19
Shaft Foam Arbors	19
Urethane and Felt Arbors	20
Pneumatic Sand Drum	20
Rubber Sand Cones	21
Tapered Felt Arbor/ Core	22
Expansion Arbors	23
Dual Purpose Expansion Arbors	24
Grinding Rasps	25
Router Adapters	25
Tycro Polishing Arbor	26
Trautman Coneiak Polishing Arbors	26

Trautman Vacuum Systems

Holding Frames	29
Drape Forming Adapters	29

PDQ Ovens

PDQ Rollabout Oven	31
PDQ BT-4 Infrared Oven	32
PDQ PO-2 Infrared Oven	32

Fabrication Equipment

Laminating Fixture LJ-100	35
Surgical Tubing	35
Vertical Fabricating Fixture, VFJ-100	36
Pipe Holding Fixtures	37
Prosthetic Bench Vise	37
Cable Swaging Tool	38
Universal Alignment Fixture	38
Contouring Fork and Bench Receptacle	39
Joint Alignment Tool	39
Bending Irons	40
Contouring Irons	40
Spring Clamps	41
Yates Clamp	41
Scissors	41
Hex Screwdriver	42
Hex Keys	42
Torque Screwdriver	42
Torque Wrench	42
Scarpas Knife	43
Surgical Scalpel	43
Pin Wrench	43
Deburring Tool	44
Cutting Tool for 30 mm Tubes	44
Cast Cutter	44
Measuring Tapes	45
Cast Mandrels	45
Surform™ Files	45
Transtibial Length Gauge	46
Original AP-ML Gauge	46
Universal AP-ML Gauge	46
Transfemoral Length Gauge	47
Rounded Blade Spatula	47

Materials

Cenpol	49
Thermo Foam 4E/ Firm Foam	50
Polypropylene	50
Thermoflex	51
Cleani-Stump Antiseptic Wipe	51
PVA Bags	52
Alginate Impression Powder	52
PolyCar-C Sheets	53

Silicone	53
Transforming Silicone Putty	54
Polyester Stockinette - Lamination Tube	55
Elastic Stockinette Tube, 7% Lycra	55
Elastic Stockinette Tube, 12,5% Lycra	56
Cotton Stockinette	56
Dacron Tape	57
Casting Socks	57
Comfort Casting Socks	57
Thermoformable Composite (TFC), Sheets	58
Thermoformable Composite (TFC), Insole	59
Bars	60
TFC Upright Bars & Silicone Forming Sleeve	60
Cotton Stockinette	61
Teflon Sheet	61
Silicone Sheet	62
Vacuum Bag material	62
Tape	62
Silicone Glove	63
Release Film - sheet	63
Silicone Bag	63

Catalog symbol legend

Specification symbols


- 

Manufactured in Sweden
- 

Manufactured in USA
- 

Manufactured in Spain


Color symbols


- 

CC03/Tan
- 

CC99/Black
- 

Grey
- 

Green
- 

Pink
- 

Orange

Trautman Carvers & Sanders

Specification symbols

SE Manufactured in Sweden

US Manufactured in USA



Trautman Double Drum Sander

US

Specifications

Height: 44 in. (112 cm)	Warranty: 12 months from the date of shipment
Width: 45 in. (114 cm)	
Weight: 535 lbs. (243 kg)	
Speed: 1725 rpm	

Product Features

The Trautman 5005 DD Sander comes with same motor, motor control, dual drive belts with guard and base as 5001 Arbor. This machine has two 9 in. (23 cm) Pneumatic Sanding Drums.

- Dual drive belts
- 2 horse power, Baldor® motor
- Durable powder coated finish

Ordering information

TRAUTMAN DOUBLE DRUM SANDER

DESCRIPTION	PRODUCT NUMBER
Trautman Arbor, Single Phase, 50Hz, 230 V	5005-1-50
Trautman Arbor, Three Phase, 50Hz, 460 V	5005-3-50

SAND CONES FOR NON-PNEUMATIC DRUMS

DESCRIPTION	PRODUCT NUMBER
Sand Drum Sleeve, 24 Grit	S24
Sand Drum Sleeve, 50 Grit	S50
Sand Drum Sleeve, 120 Grit	S120
Sand Drum Sleeve, 180 Grit	S180
Sand Drum Sleeve, 240 Grit	S240





Trautman Variable Speed Floor Carver

US

Specifications

- Height: 63 in. (160 cm)
- Working Height: 26 - 56 in. (66 - 142 cm)
- Weight: 228 lbs. (103 kg)
- Speed: 1380 - 3450 rpm, 115 V (1380 - 2860 rpm, 230 V)
- Thread: 1/2"-13
- Warranty: 12 Months from the date of shipment

Product Features

- Precise speed control for best results possible
- Heavy duty 3/4 hp Baldor® motor and controller
- Brushless motor for low maintenance
- Motor is thermally protected to prevent damage from overheating
- Factory lubricated and sealed bearings for long-lasting use
- Pre-attached 40 watt work light for the 60 Hz machine
- Fully adjustable height and tilt
- Torque overload safety feature
- Heavy duty base for stability

Ordering information

TRAUTMAN VARIABLE SPEED FLOOR CARVER

DESCRIPTION	FREQUENCY	THREAD	PRODUCT NUMBER
Trautman Floor Carver	60 Hz, 115 V	1/2"-13	2300/60
Trautman Floor Carver	50 Hz, 230 V	1/2"-13	2300/50
Anchor Bolt System			390900
Lift assist holder			2100LA
Accessory Holder			2300M



Floor Carver
shown with
Carver Lift Assist

Trautman Single Speed Floor Carver

US

Specifications

Height: 63 in. (160 cm)
Working Height: 26 - 56 in. (66 - 142 cm)
Weight: 219 lbs (99 kg)
Speed: 3450 rpm, 115 V (2860 rpm, 230 V)
Thread: 1/2"-13
Warranty: 12 Months from the date of shipment

Product Features

- Durable 1/2 hp Baldor® Motor
- Brushless motor for low maintenance
- Motor is thermally protected to prevent damage from overheating
- Factory lubricated and sealed bearings for long-lasting use
- Pre-attached 40 watt work light for the 60 Hz machine
- Fully adjustable height and tilt
- Torque overload safety feature
- 5/8"-11 threaded accessories may be used with Trautman Router Adapter, PN 225029



Floor Carver
shown with
Carver Lift Assist

Ordering information

TRAUTMAN SINGLE SPEED FLOOR CARVER

DESCRIPTION	PRODUCT NUMBER
Trautman Floor Carver, Single Speed, 50 Hz, 230 V	2100/50
Anchor Bolt System	390900
Lift Assist Holder	2100LA
Accessory Holder	2300M



Trautman Tech Station


US

Specifications

- Height: 70 in. (180 cm)
- Width: 37 in. (94 cm)
- Depth: 36 in. (91 cm)
- *Shipping Weight: 190 lbs (86 kg)
- Warranty: 12 Months from the date of shipment

Product Features

- State-of-the-art modular, multipurpose workstation
- Small footprint maximizes workspace
- 360-degree rotation with easy adjustments
- Base unit can be modified with various attachments and accessories for custom needs



Ordering information

TRAUTMAN TECH STATION

DESCRIPTION	SHIPPING WEIGHT	PRODUCT NUMBER
Trautman Tech Station	86 kg	395400

TRAUTMAN TECH STATION ATTACHMENTS

DESCRIPTION	SHIPPING WEIGHT	PRODUCT NUMBER
Trautman Tech Pipe Vise	4kg	395012
Trautman Tech Pipe Vise with Drop Hitch	6 kg	395250
Trautman Tech Drape Forming Pipe Holder	5 kg	395550
Trautman Tech Attachment for Bench Carver	3 kg	395050
Trautman 12 in. Platen	5 kg	227347
Trautman Tech LJ-100 Adapter	5 kg	395810

Dust Collector for Floor Carver

US

Product Features

- Infinitely adjustable
- Steel construction
- 4 in. (10 cm) exhaust for connection to vacuum system



Ordering information

DESCRIPTION

Dust Collector for Floor Carver

PRODUCT NUMBER

2100L

Dust Collection Hose Kit for 2100L

US

Product Features

- For use with Dust Collector for Floor Carver (PN 2100L)
- 4 in. (10 cm) diameter clear hose for viewing any obstruction
- 6 ft. (1.8 m) of flexible hose reach
- 2 worm-drive clamps for connection
- Internal abrasion resistant
- Enclosed steel wire reinforced



Ordering information

DESCRIPTION

Dust Collection Hose Kit for 2100L

PRODUCT NUMBER

226005



Air Attachment for Floor Carver

US

Product Features

The air attachment directs an adjustable air stream to blow dust and chips from the cutting tool and to assist in cooling down thermoplastic grinding for more efficient grinding. Standard 1/4 in. (0.6 cm) pipe fitting to be attached to your air supply.

- Features a flat air flow nozzle
- Flexible stay set hose for adjustable air stream



Ordering information

DESCRIPTION	PRODUCT NUMBER
Air Attachment for Floor Carver	2100J

Accessory Holder for Floor Carver

US

Product Features

- Round, rotatable carousel design
- Keeps accessories organized, visible and accessible
- 23 total peg locations, includes 8 long and 23 short pegs for custom peg layouts

NOTE: Accessories not included



Accessories sold separately

Ordering information

DESCRIPTION	PRODUCT NUMBER
Accessory Holder for Floor Carver	2300M

Lift Assist Kit for Floor Carver

US

Product Features

- Assists in both lifting and lowering
- Easy to install
- Provides over 16 in. (40 cm) of adjustment
- Works on all Trautman models, past and current
- Thrust bearing for ease of rotation
- Shock provides slam-resistant damping at the end of the stroke



Ordering information

DESCRIPTION	PRODUCT NUMBER
Carver Lift Assist	2100LA



Trautman Variable Speed Bench Carver

US

Specifications

Height: 22 in. (56 cm)	Speed: 1380 - 3400 rpm, 115 V (1380 - 2820 rpm, 230 V)
Working Height: 15 - 22 in. (38 - 56 cm)	Thread: 1/2"-13
Weight: 28 lbs. (12.7 kg)	Warranty: 12 months from the date of shipment

Product Features

- Maintenance-free 1/3 hp motor and controller
- High quality variable speed motor controller
- Lubricated ball bearings, sealed for life

Ordering information

TRAUTMAN ARBOR

DESCRIPTION	PRODUCT NUMBER
Trautman Bench Carver, Variable Speed, 50 Hz, 230 V, 1/2"-13 Thread	2600/50



Trautman Single Speed Bench Carver

US

Specifications

Height: 22 in. (56 cm)	Speed: 3400 rpm, 115 V (2820 rpm, 230 V)
Working Height: 15 - 22 in. (38 - 56 cm)	Thread: 1/2"-13
Weight: 25 lbs. (11.3 kg)	Warranty: 12 Months from the date of shipment

Product Features

- Great for satellite offices or bench work
- Maintenance free 1/3 hp motor
- Lubricated sealed ball bearings
- 5/8-11 threaded accessories may be used with Trautman Router Adapter, PN 225029)

Ordering information

TRAUTMAN ARBOR

DESCRIPTION	PRODUCT NUMBER
Trautman Bench Carver, 50 Hz, 230 V, 1/2"-13 Thread	2500/50



Bench Carver Floor Stand

US

Product Features

- Adjustable height 35 – 56 in. (89 - 142 cm)
- Alloy screw adjusts clamping force



Ordering information

DESCRIPTION

Bench Carver Floor Stand

PRODUCT NUMBER

390190

Accessory Holder for Bench Carver

US

Product Features

- Keeps accessories organized, visible, and accessible



Ordering information

DESCRIPTION

Accessory Holder for Bench Carver

PRODUCT NUMBER

2500M

Dust Collector for Bench Carver

US

Product Features

- Works with most workshop vacuum systems
- Flexible gooseneck for wide range of adjustment
- Appropriate for bench carvers and for hand-held sanding tools



Ordering information

DESCRIPTION

Dust Collector for Bench Carver

PRODUCT NUMBER

2500L



Drill Attachments 1/2"-13 Thread

US



Ordering information

PRODUCT FEATURES	PRODUCT NUMBER
Trautman Drill Attachment - Fits up to 3/8 in. shaft	2100S

Soft Foam Shaper

US

Specifications

Length: 5 in. (12 cm)

Thread: 5/8"-11

Max. Diameter Schmid Carver: 2-3/4 in.

Product Features

- Ideal for contouring both inner and outer areas of foam cover
- Fine steel bristles barely protrude through surface so foam material may be shaped quickly and removed to obtain the desired cosmetic finish
- Trautman Router Adapter required (PN 225029)



Ordering information

DESCRIPTION	PRODUCT NUMBER
Soft Foam Shaper	51391

Shaft Foam Arbors

US

Ordering information

FOAM ARBOR FOR FLOOR CARVER

DESCRIPTION	PRODUCT NUMBER
Shaft Foam Arbor for Floor Carver, 3 × 9 in.	2100X

FOAM ARBOR FOR BENCH CARVER

DESCRIPTION	PRODUCT NUMBER
Shaft Foam Arbor for Bench Carver, 3 × 5-1/2 in.	2500X

2100X SAND CONES

DESCRIPTION	PRODUCT NUMBER
Sand Cone for 2100X, 36 Grit	225240
Sand Cone for 2100X, 60 Grit	225242
Sand Cone for 2100X, 120 Grit	225244

2500X SAND CONES

DESCRIPTION	PRODUCT NUMBER
Sand Cone for 2500X, 36 Grit	225230
Sand Cone for 2500X, 60 Grit	225232
Sand Cone for 2500X, 120 Grit	225234



Urethane and Felt Arbors

US

Product Features

- **Urethane Arbors:** Most commonly used for shaping soft foam covers
- **Felt Arbors:** Used for polishing plastic edges
- Can be shaped to any contour



225080

Ordering information

URETHANE FEATURES

DESCRIPTION	PRODUCT NUMBER
Urethane Arbor, 3 × 5 in., 5/8"-11 Thread	225080
Urethane Arbor, 3 × 5 in., 1/2"-13 Thread	225082

Pneumatic Sand Drum

US

Ordering information

PNEUMATIC SAND DRUM 3.5 X 5 IN.

DESCRIPTION	PRODUCT NUMBER
Pneumatic Sand Drum, 1/2"-13 Thread	224100



SAND DRUM SLEEVES - PNEUMATIC

DESCRIPTION	PRODUCT NUMBER
Sand Sleeve, 120 Grit	224101
Sand Sleeve, 80 Grit	224102
Sand Sleeve, 36 Grit	224103

Rubber Sand Cones

US

Specifications

Length: 4-3/8 in. (111 mm)

Thread: 1/2"-13

Ordering information

Mini Sand Core



MINI SAND CORE

DESCRIPTION	PRODUCT NUMBER
Sand Core, Mini	2100P

MINI SANDING SLEEVES

DESCRIPTION	PRODUCT NUMBER
Mini 50 Grit	P50
Mini 100 Grit	P100
Mini 180 Grit	P180

MINI POLISHING CONE

DESCRIPTION	PRODUCT NUMBER
Mini, Fine, Blue	225117
Mini, Medium, Maroon	225118
Mini, Coarse, Brown	225119

Small Sand Core



SMALL SAND CORE

DESCRIPTION	PRODUCT NUMBER
Sand Core, Small	2100D

SMALL SANDING SLEEVES

DESCRIPTION	PRODUCT NUMBER
Small 24 Grit	D24
Small 80 Grit	D80
Small 120 Grit	D120

SMALL POLISHING CONE

DESCRIPTION	PRODUCT NUMBER
Small, Fine, Blue	225111
Small, Medium, Maroon	225112

Medium Sand Core



MEDIUM SAND CORE

DESCRIPTION	PRODUCT NUMBER
Sand Core, Medium	2100C

MEDIUM SANDING SLEEVES

DESCRIPTION	PRODUCT NUMBER
Medium 24 Grit	C24
Medium 80 Grit	C80
Medium 120 Grit	C120
Medium 180 Grit	C180
Medium 240 Grit	C240
Medium 320 Grit	C320

MEDIUM POLISHING CONE

DESCRIPTION	PRODUCT NUMBER
Medium, Fine, Blue	225114
Medium, Medium, Maroon	225115
Medium, Coarse, Brown	225116



Tapered Felt Arbor/ Core

US



Ordering information

1×4 IN. FELT ARBORS & SAND CONES

DESCRIPTION	PRODUCT NUMBER
1 × 4 in. Sand Cone, 80 Grit	225173
1 × 4 in. Sand Cone, 120 Grit	225175

1.5×4 IN. FELT ARBORS & SAND CONES

DESCRIPTION	PRODUCT NUMBER
1.5 × 4 in. Sand Cone, 36 Grit	225151
1.5 × 4 in. Sand Cone, 100 Grit	225169
1.5 × 4 in. Sand Cone, 120 Grit	225171

2×4 IN. FELT ARBORS & SAND CONES

DESCRIPTION	PRODUCT NUMBER
2 × 4 in. Sand Cone, 36 Grit	225136
2 × 4 in. Sand Cone, 80 Grit	225141



Expansion Arbors

US

Specifications

Warranty: 6 months against manufacturing defects

Product Features

- Expands to apply inside pressure to the sanding sleeve to hold it in place
- Prevents migration and “Fly away” of sanding sleeve
- Provides a smoother grinding surface for straighter lines and higher quality finish
- Eliminates grinding chatter that hard grinding surfaces gives
- Sold individually and as a kit of three
- Available in 1/2"-13 and M16 threads to fit most O&P machinery



Ordering information

EXPANSION ARBOR KITS

DESCRIPTION	PRODUCT NUMBER
1/2"-13 Threads	225380
M16 Threads	225380MM

Mini



MINI EXPANSION ARBOR

DESCRIPTION	PRODUCT NUMBER
1/2"-13 Thread	225350
M16 Thread	225350MM

MINI SANDING SLEEVES

DESCRIPTION	PRODUCT NUMBER
Mini 50 Grit	P50
Mini 100 Grit	P100
Mini 180 Grit	P180

Small



SMALL EXPANSION ARBOR

DESCRIPTION	PRODUCT NUMBER
1/2"-13 Thread	225352
M16 Thread	225352MM

SMALL SANDING SLEEVES

DESCRIPTION	PRODUCT NUMBER
Small 24 Grit	D24
Small 80 Grit	D80
Small 120 Grit	D120

Medium



MEDIUM EXPANSION ARBOR

DESCRIPTION	PRODUCT NUMBER
1/2"-13 Thread	225354
M16 Thread	225354MM

MEDIUM SANDING SLEEVES

DESCRIPTION	PRODUCT NUMBER
Medium 24 Grit	C24
Medium 80 Grit	C80
Medium 120 Grit	C120
Medium 180 Grit	C180 (metal use)
Medium 240 Grit	C240 (metal use)
Medium 320 Grit	C320 (metal use)



Dual Purpose Expansion Arbors

US

Specifications

Warranty: 6 months against manufacturing defects

Product Features

- Sand and polish without changing tools
- Removable, tapered polishers for easy replacement

Ordering information

Small



SMALL DUAL PURPOSE

DESCRIPTION	PRODUCT NUMBER
1/2"-13 Threads	225330
5/8"-11 Threads	225335
M16 Threads	225330MM
Slim Blue Coneiak, Replacement tip	225345

SMALL SANDING SLEEVES

DESCRIPTION	PRODUCT NUMBER
Small 24 Grit	D24
Small 80 Grit	D80
Small 120 Grit	D120

Medium



MEDIUM DUAL PURPOSE

DESCRIPTION	PRODUCT NUMBER
1/2"-13 Threads	225320
5/8"-11 Threads	225325
M16 Threads	225320MM
Tapered Red Coneiak, Replacement tip	225340

MEDIUM SANDING SLEEVES

DESCRIPTION	PRODUCT NUMBER
Medium 24 Grit	C24
Medium 80 Grit	C80
Medium 120 Grit	C120
Medium 180 Grit	C180 (metal use)
Medium 240 Grit	C240 (metal use)
Medium 320 Grit	C320 (metal use)

Grinding Rasps

US

Product Features

Warranty: 6 months against manufacturing defects

- High-quality steel
- Sharp teeth
- Long usage

Ordering information

5/8"-11 THREAD*

DESCRIPTION	PRODUCT NUMBER
1-1/2 × 2-1/4 in. - Use on Transtibial Sockets, Coarse	225052
1-3/8 × 3-3/8 in. - Small Diameter and Channel Work, Coarse	225045

*Must have router adapter for 5/8-11 threaded Rasps


1/2"-13 THREAD

DESCRIPTION	PRODUCT NUMBER
Rasp	2100E
Carbide Grind Rasp, Large, Fine	225256



225052

225045



2100E

225256

Router Adapters

US

Specifications

Warranty: 6 months against manufacturing defects

Ordering information

1/2"-13 INTERNAL THREAD

DESCRIPTION	PRODUCT NUMBER
5/8-11 External Thread	225029
M16 External Thread	225017

1/2"-13 EXTERNAL THREAD

DESCRIPTION	PRODUCT NUMBER
5/8-11 Internal Thread	225011



225017

225029



225011

Tycro Polishing Arbor

US

Specifications

Warranty: 6 months against manufacturing defects



225100

Ordering information

TYCRO ACCESSORIES

DESCRIPTION	PRODUCT NUMBER
Tycro Cone, 4 x 1 in. (Requires 225090)	225100
Tycro Cone Adapter, 1/2 in. Thread	225090

Trautman Coneiak Polishing Arbors

US

Specifications

Warranty: 6 months against manufacturing defects



225600 225615 225605 225620 225610

Product Features

- All Coneiak Polishing Arbors work well with any standard thermoplastic materials within the O&P industry.
- For best results, drill attachment 2100S (large) or 2100R (small) should be used
- If using a rotary hand tool, the maximum speed should not exceed 10,000 RPMs
- Ensure that the shaft is at least halfway into chuck

Ordering information - Coneiak Mini Polishing Arbor

DESCRIPTION	SHAFT	PRODUCT NUMBER
Set of Five Coneiak Minis	1/8 in.	225650
Tapered Coneiak Mini	1/8 in.	225600
Ball Nose Coneiak Mini	1/8 in.	225615
Slim Coneiak Mini	1/8 in.	225605
Cylindrical Coneiak Mini	1/8 in.	225620
Double Ended Coneiak Mini	1/8 in.	225610

Ordering information - Trautman Coneiak Polishing Arbor

US

SET OF FOUR CONEIAKS

DESCRIPTION	PRODUCT NUMBER
1/2"-13 Thread	225490
5/8"-11 Thread	225495
M16 Thread	225490MM



TAPERED CONEIAK

DESCRIPTION	PRODUCT NUMBER
1/2"-13 Thread	225410
5/8"-11 Thread	225415
M16 Thread	225410MM



BALL NOSE CONEIAK

DESCRIPTION	PRODUCT NUMBER
1/2"-13 Thread	225400
5/8"-11 Thread	225405
M16 Thread	225400MM



SLIM CONEIAK

DESCRIPTION	PRODUCT NUMBER
1/2"-13 Thread	225440
5/8"-11 Thread	225445
M16 Thread	225440MM



CYLINDRICAL CONEIAK

DESCRIPTION	PRODUCT NUMBER
1/2"-13 Thread	225450
5/8"-11 Thread	225455
M16 Thread	225450MM



HOURLASS CONEIAK

DESCRIPTION	PRODUCT NUMBER
1/2"-13 Thread	225430
5/8"-11 Thread	225435
M16 Thread	225430MM



GROOVED CONEIAK

DESCRIPTION	PRODUCT NUMBER
1/2"-13 Thread	225420
5/8"-11 Thread	225425
M16 Thread	225420MM



DOUBLE DUROMETER

DESCRIPTION	PRODUCT NUMBER
1/2"-13 Thread	225460
5/8"-11 Thread	225465
M16 Thread	225460MM





Trautman Vacuum Systems

Specification symbols

SE Manufactured in Sweden

US Manufactured in USA

Holding Frames

US


Product Features

- Black powder coated for durability

Ordering information

Note: Sold per pair

DESCRIPTION	PRODUCT NUMBER
Holding Frames, 12 × 12 in.	227420
Holding Frames, 18 × 18 in.	227443
Holding Frames, 20 × 20 in.	227449
Holding Frames, 24 × 24 in.	227451



Drape Forming Adapters

US



Ordering information

PLATEN

DESCRIPTION	PRODUCT NUMBER
9 in. Platen	227354
12 in. Platen	227347

PDQ Ovens

Specification symbols

SE Manufactured in Sweden

US Manufactured in USA

PDQ Rollabout Oven

US

Specifications

Electrical requirements: US 208 - 240 V single phase
EU 380 - 460 V three phase

Warranty: 12 months from the date of shipment

Product Features

A Rollabout model is an excellent solution for heating large sheets of plastic for hip. There are 4 different sizes to meet every thermoforming need, and each features a hydraulic assisted height adjustable tray for easy operation. The tray includes a fitted Teflon® cover, is mounted on casters, and can be rolled to your vacuum area. No more burnt arms or trying to find a place to set down your hot plastic.

- Aluminum and stainless steel construction
- 650 watt infrared emitters (1000 watt on 8E & 8D)
- Double pane viewing window & interior light
- Programmable digital temperature controller
- Tray includes a fitted Teflon® cover
- NOTE: 208 - 240 V, single phase ovens delivered without cord
380 - 460 V, three phase ovens delivered with cord, no plug

Dimensions

MODEL	INTERIOR DIMENSIONS (W X D X H)	EXTERIOR DIMENSIONS (W X D X H)	WEIGHT
PDQ-5A	38 × 31 × 18 in. (965 x 787 x 457 mm)	49 × 35 × 48 in. (1244 x 889 x 1219 mm)	350 lbs. (159 kg)
PDQ-6B	36 × 51 × 18 in. (914 x 1295 x 457 mm)	47 × 55 × 48 in. (1194 x 1397 x 1219 mm)	450 lbs. (204 kg)
PDQ-7C	48 × 51 × 18 in. (1219 x 1295 x 457 mm)	59 × 55 × 48 in. (1499 x 1397 x 1219 mm)	460 lbs. (209 kg)
PDQ-8E	60 × 50 × 18 in. (1524 x 1219 x 457 mm)	70 × 52 × 48 in. (1778 x 1321 x 1219 mm)	475 lbs. (215 kg)

Ordering information

MODEL	PROD NO. EU 400 V, 50 HZ, THREE PHASE
PDQ-5A	O5A-INT
PDQ-6B	O6B-INT
PDQ-7C	O7C-INT
PDQ-8E	O8E-INT



Shown with the PDQ Oven Fan Kit.
See page 37 for more information.



PDQ BT-4 Infrared Oven

US

Specifications

Electrical requirements: US 208 - 240 V single phase
EU 380 - 460 V three phase

Warranty: 12 months from the date of shipment

Product Features

The BT-4 is capable of bubble forming and heating sheet plastic while taking up a minimal amount of floor space. The BT-4 Infrared Oven is big enough to heat most KAFO sheet plastics and has enough drop for even the biggest check socket.

The BT-4 has the same speed and reliability you’ve come to expect from PDQ Infrared Ovens. The eight 1000 watt emitters, four on top and four on bottom, along with the dual electric fans allow bubble forming plastics to be heated evenly and efficiently. The BT-4 is truly the best of both worlds with its small size and large capacity.

- Eight 1000 watt emitters; four on top and four on bottom
- Programmable digital temperature controller
- Dual electric fans for air circulation
- Access point for vacuum line
- 4 interior lights and large viewing window
- Sheet plastic tray with fitted Teflon® cover
- Aluminum and powder coated protective finish
- 18 in. (46 cm) clearance for bubble forming
- NOTE: US 220 V, single phase ovens delivered without cord
EU 440 V, three phase ovens delivered with cord, no plug



Ordering information

MODEL	INTERIOR DIMENSIONS (W x D x H)	EXTERIOR DIMENSIONS	PRODUCT NUMBER
PDQ-BT-4	27 × 40 × 24 in. (686 x 1016 x 610 mm)	32 × 49 × 61 in. (813 x 1245 x 1549 mm)	o-BT4-INT

PDQ PO-2 Infrared Oven

US

Specifications

Electrical requirements: US 208 - 240 V single phase
EU 380 - 460 V three phase

Warranty: 12 months from the date of shipment

Product Features

Designed with banks of lights on top and bottom to facilitate the quick and even forming of thick plastics.

Instead of having to reach into a hot oven for bubble forming frames, the PO-2 model now has a rollabout tray that holds either frames or trays for heating flat plastic. Just roll the cart out to your mold.

- Aluminum and stainless steel construction
- Quartz infrared tubes top and bottom
- Eight 1000 watt quartz infrared tubes (four on top and four on bottom)
- Interior light
- Two double pane viewing windows
- Programmable digital temperature controller (accurate to within 2°)
- Rollabout bubble forming rack/sheet plastic tray with fitted Teflon® cover
- Bubble Forming Frame (specify size at time of ordering)
- Full 1 year warranty
- NOTE: US 220 V, singel phase ovens delivered without cord
EU 440 V, three phase ovens delivered with cord, no plug

Ordering information

MODEL	INTERIOR DIMENSIONS (W x D x H)	EXTERIOR DIMENSIONS	PRODUCT NUMBER
PDQ-PO-2	50 × 36 × 20 in. (1270 x 914 x 508 mm)	60 × 40 × 53 in. (1524 x 1016 x 1346 mm)	O-PO2-INT





Fabrication Equipment

Specification symbols

- SE Manufactured in Sweden
- US Manufactured in USA



Laminating Fixture LJ-100

US

Product Features


- Height: 24 in. (61 cm)
- Width: 17 in. (43 cm)
- May be positioned to hold cast at any angle or degree of rotation
- Large PVA tie-off ring and hose attachments for inside and outside vacuum lines make laminating easier and faster

Ordering information

DESCRIPTION	PRODUCT NUMBER
Laminating Fixture Complete, LJ-100	51454
Vacuum Device Assembly Only	51459
Vacuum Laminating Assembly Only	51455

SPARE PARTS FOR LJ100

A.	51462 - Clamp Assembly
B.	51463 - Post Base Assembly
C.	52062 - Clamp Assembly
D.	51271 - Hose Connector
E.	51458 - Stopper
F.	51457 - Plug
G.	51456 - Pipe Lock
H.	Surgical Tubing (See Below)
I.	55706 - Lamination Mandrell
	51275 - Hose, Tie-Off 3/16 in., (Not Shown)



51459 Vacuum Laminating Device

Surgical Tubing

US

Product Features

- High quality, pure latex tubing
- May be used as hose on vacuum or to securely tie-off PVA bags on the Vacuum Laminating Fixture, LJ-100

Ordering information

DESCRIPTION	PRODUCT NUMBER
Surgical Tubing, 0.25 in. (6.4 mm) Diameter, 6 ft. (183 cm) Increments Only	51273
Surgical Tubing, 0.19 in. (4.8 mm) Diameter, per foot	51275





Vertical Fabricating Fixture, VFJ-100

US

Product Features

- **Height:** 48-3/4 in. (124 cm)
- **Width:** 21 in. (53 cm)
- Horizontal operation
- Heavy-duty for long life
- Retractable handle when not in use
- Will form full circle without restriction



Ordering information

DESCRIPTION	PRODUCT NUMBER
Vertical Fabricating Fixture, VFJ-100	52056
Vertical Swivel Module	61714
VFJ-Ball Joint Module	58188

Pipe Holding Fixtures

US

Ordering information

DESCRIPTION	PRODUCT NUMBER
Pipe Holding Fixture, Vice Mount	220129
Bench Mounted Pipe Holding Fixture, 8 in. Hi	220131
Friction Block	220135
Laminating Stand with Pipe Holding Fixture Adjustable in Height to 60 in. (152 cm), 15 – 18 in. (381 - 457 mm) Diameter	220133



Prosthetic Bench Vise

US

Product Features

- Tough, resilient pads
- Opens to 16 in. (406 mm) with an 8 in. (203 mm) depth
- Swivel jaw, 6 × 6 in. (152 x 152 mm)

Ordering information

DESCRIPTION	PRODUCT NUMBER
Prosthetic Bench Vise	220178
Replacement Pads	220301





Cable Swaging Tool

US

Product Features

- May be used in any standard bench vise
- Made of heat treated steel

Ordering information

DESCRIPTION	PRODUCT NUMBER
Cable Swaging Tool, for Three Sizes of Cable Size Fittings	50520

REPLACEMENT DIE

DESCRIPTION	PRODUCT NUMBER
Replacement Die, 3/64 in. (0.119cm)	50524
Replacement Die, 1/16 in. (0.159cm)	50523
Replacement Die, 3/32 in. (0.238cm)	50522



Universal Alignment Fixture

US

Ordering information

DESCRIPTION	PRODUCT NUMBER
Universal Alignment Fixture	220128



Contouring Fork and Bench Receptacle


US

Product Features

- Heat-treated stainless steel
- Mounts flush on work bench or in work bench vise

Ordering information

DESCRIPTION	PRODUCT NUMBER
Contouring Fork, 5/16 in. Throat	220194



Joint Alignment Tool


US

Product Features

- Use for both small and large stirrups

Ordering information

DESCRIPTION	PRODUCT NUMBER
Alignment Tool, Universal 0.25 in. diameter × 8 in. long rod	220194





Bending Irons

US

Product Features

- 8 in. (203 mm) long, 1/4 in. (6.4 mm) thickness capacity
- Stainless steel
- Provides smooth bending of side bars
- Sold as a pair



Ordering information

DESCRIPTION	PRODUCT NUMBER
Bending Irons	54795

Contouring Irons

US

Product Features

- Nickel plated bronze
- Wide-radius jaws
- Sold as a pair
- Will not mark finished metal or leather work



Ordering information

DESCRIPTION	PRODUCT NUMBER
Rounded Polished Jaws	220285

Spring Clamps

US



Ordering information

DESCRIPTION	PRODUCT NUMBER
Spring Clamps Set (4 Each)	227594

Yates Clamp

US

Product Features

- Holds harness during fitting



Ordering information

DESCRIPTION	LENGTH	PRODUCT NUMBER
Wire Spring Type, W1 Small	3.1 in. (8 cm)	52115
Wire Spring Type, W2 Medium	3.7 in. (9.5 cm)	52116
Rubber Pads, for 52115 and 52116		52872

Scissors



5-514

Ordering information

DESCRIPTION	PRODUCT NUMBER
Bandage Scissors, 5.5 in. (140 mm)	5-514
Bandage Scissors, 7.25 in. (184 mm)	5-516
Bandage Scissors, 8 in. (203 mm)	5-550



Tools

US

Hex Screwdriver




DESCRIPTION	PRODUCT NUMBER
Hex Screwdriver, for U-30/G-2 Adjustment Screws, 4 mm Hex	731000120
Hex Screwdriver, for U-30/G-2 Fixering Screws, 3 mm Hex	731000121

Hex Keys

US

SE



DESCRIPTION	PRODUCT NUMBER
Hex Key, 4 mm	51293
Hex Key, 3 mm	51353

Torque Screwdriver

SE



DESCRIPTION	PRODUCT NUMBER
BOA Torque Screwdriver, 1.6 Nm	731000153
1/4 in. Bits 4 mm Hex, for BOA Torque Screwdriver	731000155

Torque Wrench

SE



DESCRIPTION	PRODUCT NUMBER
BOA Torque Wrench, 15 Nm	731000150

Scarpas Knife

US



Ordering information

DESCRIPTION	PRODUCT NUMBER
Scarpas Knife, Small for Transtibial	220461
Scarpas Knife, Large for Transfemoral	220460

Surgical Scalpel

US



Ordering information

DESCRIPTION	PRODUCT NUMBER
Scalpel Handle	800050

Pin Wrench

Product Features

- Used for tightening: ALP-Ti Female Pyramid Adapter, 731162120 and ALP-Al Female Pyramid Adapter, 731163187



Ordering information

DESCRIPTION	PRODUCT NUMBER
Pin Wrench DIN1810B-45-50	731000140



Deburring Tool

US

Ordering information



DESCRIPTION	PRODUCT NUMBER
Deburring Tool	220321
Deburring Tool, Blade Only	E-100

Cutting Tool for 30 mm Tubes

US

Ordering information



DESCRIPTION	PRODUCT NUMBER
Cutting Tool, for 30 mm Carbon Fiber Tube	731000125

Cast Cutter

US

Product Features

- Manufacturer’s descriptive literature available on request
- Oscillating action quickly cuts through plaster
- Trim plastic sockets and Ortholen splints

Ordering information



DESCRIPTION	PRODUCT NUMBER
Cast Cutter Blade, 2.5 in. (64 mm)	226258

Measuring Tapes

US



220331

220328

220329

Ordering information

DESCRIPTION	PRODUCT NUMBER
Measuring Tape (Inch/Metric) Fiberglass Retractable, 60 in. (150 cm)	220328
Measuring Tape (Inch/Metric) PVC-Coated, Fiberglass 60 in. (150 cm)	220329
Disposable Measuring Tape, 10 ft. (300 cm)	220331

Cast Mandrels

US



Ordering information

DESCRIPTION	PRODUCT NUMBER
Aluminum Cast Mandrel, 1.75 × 7.5 in. (44 x 190 mm)	259056
Aluminum Cast Mandrel, 2.38 × 7.5 in. (61 x 190 mm)	259060
Aluminum Cast Mandrel, 2.38 × 10 in. (61 x 254 mm)	259053
Aluminum Cast Mandrel, 2.38 × 15 in. (61 x 381 mm)	259050

Surform™ Files

US

Product Features

- For fast shaping of plaster cast
- Also can be used on wood and plastics

Ordering information

DESCRIPTION	PRODUCT NUMBER
Flat File 10 in. (254 mm), Complete	21-295

REPLACEMENT BLADES

DESCRIPTION	PRODUCT NUMBER
Flat Blade	21-293
Half Round Blade	21-299





Transtibial Length Gauge

US

Product Features

- Measures in inches and metric



Ordering information

DESCRIPTION	PRODUCT NUMBER
Transtibial Gauge	220145

Original AP-ML Gauge

US

Product Features

- Anatomical measuring tool
- Measures in inches and metric



Ordering information


DESCRIPTION	PRODUCT NUMBER
Original AP-ML Gauge, 22 in.	220335

Universal AP-ML Gauge

US

Product Features

- Anatomical measuring tool
- Measures in inches and metric



Ordering information

DESCRIPTION	PRODUCT NUMBER
Universal Gauge, with Level, 17 in.	220336

Transfemoral Length Gauge

US

Product Features

- Measures in inches and metric



Ordering information

DESCRIPTION	PRODUCT NUMBER
Transfemoral Length Gauge	220137

Rounded Blade Spatula

US

Ordering information



DESCRIPTION	PRODUCT NUMBER
Rounded Blade Spatula 0.75 × 4 in. (19 x 102 mm) ,small	220059
Rounded Blade Spatula 1 × 6 in. (25 x 152 mm), large	220063



Materials

Specification symbols

SE Manufactured in Sweden

US Manufactured in USA

Cenpol

SE

Specifications

Density: 30 kg/cubic meter

Shore: 00=44

Product Features

- Suitable for orthoses padding and the build-up of BK Cosmesis
- Made of EVA
- Closed cell material that prevents moisture absorption
- Color: beige
- Forming temperature: 266 - 320 °F (130 - 160° C)

Ordering information

DESCRIPTION	TIME IN OVEN	PRODUCT NUMBER
39 x 39 x 0.78 in. (1000 x 1000 x 20 mm)	6 min.	637100020
39 x 39 x 1.0 in. (1000 x 1000 x 25 mm)	8 min.	637100025



Thermo Foam 4E/ Firm Foam

US
 SE

Product Features

- Closed cell material that prevents moisture absorbtion
- Thermo foam: hardness shore 00 - 58
- Firm Foam: hardness shore 00 - 88 is considerably harder and suitable for the tongue on braces
- Forming temperature 300°- 356°F (150° - 180°C)
- Color: White

Ordering information

THERMOFOAM

DESCRIPTION	PRODUCT NUMBER
866 x 27.6 x 0.12 in. (22000 x 700 x3 mm)	637700003
591 x 27.6 x 0.24 in. (15000 x 700 x 6 mm)	637700006
36 x 30 x 0.47 in. (920 x 770 x 12 mm)	637770012

FIRMFOAM

DESCRIPTION	PRODUCT NUMBER
295 x 29.9 x 0.12 in. (7500 x 760 x 3 mm)	637760003
47 x 30 x 0.47 in. (1200 x 770 x 12 mm)	637770112

Polypropylene

SE


Product Features

- Copolymer
- Durable, easy to form thermoplastic
- Forming temperature 356°- 392°F (180° - 200°C)
- Color: Natural

Ordering information

COPOLYMER POLYPROPYLENE

DESCRIPTION	PRODUCT NUMBER
48 x 32 x 0.08 in. (1220 x 800 x 2 mm)	639122002
48 x 32 x 0.12 in. (1220 x 800 x 3 mm)	639122003
48 x 32 x 0.18 in. (1220 x 800 x 4.5 mm)	639122045
48 x 32 x 0.23 in. (1220 x 800 x 6 mm)	639122006



Thermoflex

SE

Product Features

- For flexible sockets and similar
- Skin friendly transparent thermoplastic that does not shrink
- Durable material that is flexible and easy to form to provide good attachment to the skin and seating comfort
- Forming temperature 300°- 356°F (150° - 180°C)



Ordering information

DESCRIPTION	PRODUCT NUMBER
16 x 16 x 0.31 in. (400 x 400 x 8 mm), Black	637400308
16 x 16 x 0.35 in. (400 x 400 x 9 mm)	637400009
16 x 16 x 0.39 in. (400 x 400 x 10 mm), Black	637400310
16 x 16 x 0.47 in. (400 x 400 x 12 mm)	637400012
16 x 16 x 0.47 in. (400 x 400 x 12 mm), Black	637400312

Cleani-Stump Antiseptic Wipe

SE

Product Features

- Developed for use on the sensitive skin of residual limbs and can also be used to disinfect the inside of sockets
- Contains an effective alcohol and lactic acid antiseptic which reduce the amount of bacteria and some viruses on the skin for up to 16 hours
- In addition to the anti-bacterial ingredients found in Cleani-Stump, a special emollient has been added to smoothen and protect skin, helping to prevent breakdown in sensitive areas



Ordering information

DESCRIPTION	PRODUCT NUMBER
Cleani-Stump Antiseptic Wipe case (case of 40 consumer boxes, 2000 total pieces)	631020000
Cleani-Stump Antiseptic Wipe (1 consumer box, 50 wipes pieces)	631020001



PVA Bags

SE

Product Features

- Prefabricated PVA bags made of polyvinyl alcohol film
- Manufactured with a narrow, consistent seam
- 20 pieces per box



Ordering information

PVA BAGS

DESCRIPTION	PRODUCT NUMBER
4 in. (102 x 10 x 5 cm)	631100004
6 in. (102 x 15 x 5 cm)	631100006
8 in. (102 x 20 x 5 cm)	631100008
10 in. (102 x 25 x 5 cm)	631100010
12 in. (102 x 30 x 5 cm)	631100012

Alginate Impression Powder

SE

Product Features

- Ideal for casting hands, feet or reproducing sockets
- Easily dissolved in water



Ordering information


DESCRIPTION	PRODUCT NUMBER
Box 13.2 lbs. (6 kg)	633110610
Plastic Bucket 4.4 lbs. (2 kg)	633110020

PolyCar-C Sheets

US

Product Features

- High performance reinforcing composite
- Easy to mold
- High carbon content
- Thermobonds to polypropylene, copolymer and metal



Ordering information

DESCRIPTION	PRODUCT NUMBER
12 X 12 X 0.06 In. (300 X 300 X 1.5 Mm), Bidirectional Fiber Orientation	076608
12 X 24 X 0.06 In. (300 X 600 X 1.5 Mm), Bidirectional Fiber Orientation	076610
12 X 12 X 0.09 In. (300 X 300 X 2.3 Mm), Bidirectional Fiber Orientation	076612
12 X 24 X 0.09 In. (300 X 600 X 2.3 Mm), Bidirectional Fiber Orientation	076614
12 X 12 X 0.125 In. (300 X 300 X 3 Mm), Bidirectional Fiber Orientation	076616
12 X 24 X 0.125 In. (300 X 600 X 3 Mm), Bidirectional Fiber Orientation	076618
12 X 18 X 0.125 In. (300 X 450 X 3 Mm), Unidirectional Fiber Orientation	076620
12 X 36 X 0.125 In. (300 X 900 X 3 Mm), Unidirectional Fiber Orientation	076622

Silicone

SE

Ordering information

DESCRIPTION	PRODUCT NUMBER
Silicone Adhesive 90 MI	632040090





Transforming Silicone Putty



Product Features

- 1 kg - 500 g component A and 500 g component B
- Mixing ratio: 1:1
- Hardness shore A - 35



Ordering information

DESCRIPTION	PRODUCT NUMBER
1 kg Transforming Silicone Putty (500 g component A and 500 g component B)	632060010

Polyester Stockinette - Lamination Tube

SE

Product Features

- Made of 100% polyester
- Very strong and flexible, does not need to be combined with other materials to reach sufficient strength
- 22 yd. (20 m) per role



Ordering information

DESCRIPTION	WIDTH	PRODUCT NUMBER
Lamination Tube	1.2 in. (30 mm)	636203003
Lamination Tube	2.4 in. (60 mm)	636206005
Lamination Tube	3.1 in. (80 mm)	636208005
Lamination Tube	3.9 in. (100 mm)	636210010

Elastic Stockinette Tube, 7% Lycra

SE

Product Features

- Made of 93% cotton and 7% lycra
- Excellent for casting patients
- 22 yd. (20 m) per role



Ordering information

DESCRIPTION	WIDTH	PRODUCT NUMBER
Elastic Tube	1.6 in. (40 mm)	636304020
Elastic Tube	2.4 in. (60 mm)	636306020
Elastic Tube	3.1 in. (80 mm)	636308020
Elastic Tube	4.3 in. (110 mm)	636311020



Elastic Stockinette Tube, 12,5% Lycra

SE

Product Features

- Made of 87,5% cotton and 12,5% lycra
- Soft, with powerful compression
- Excellent for use under orthoses
- 22 yd. (20 m) per role



Ordering information

DESCRIPTION	WIDTH	PRODUCT NUMBER
Elastic Stockinette	1.6 in. (40 mm)	636404020
Elastic Stockinette	2.4 in. (60 mm)	636406020
Elastic Stockinette	3.1 in. (80 mm)	636408020
Elastic Stockinette	4.3 in. (110 mm)	636411020

Cotton Stockinette

SE

Product Features

- Made of 100% cotton
- Excellent for use under orthosis, casting patients, or for donning prosthesis
- 22 yd. (20 m) per role



Ordering information

DESCRIPTION	WIDTH	PRODUCT NUMBER
Cotton Stockinette	1.6 in. (40 mm)	636504020
Cotton Stockinette	2.4 in. (60 mm)	636506020
Cotton Stockinette	3.1 in. (80 mm)	636508020
Cotton Stockinette	4.3 in. (110 mm)	636511020

Dacron Tape

US

Product Features

- Durable tape has “web” type of construction to prevent stretching
- 25 yd. (22.5 m) per roll
- Available in white, black and beige



Ordering information

DESCRIPTION	WIDTH	PRODUCT NUMBER - WHITE	PRODUCT NUMBER - BEIGE	PRODUCT NUMBER - BLACK
Dacron Tape	0.5 in. (13 mm)	212167	212242	212243
Dacron Tape	0.75 in. (19 mm)	212183	212252	-
Dacron Tape	1 in. (25 mm)	212209	212245	212246
Dacron Tape	1.5 in. (38 mm)	212225	212247	212248
Dacron Tape	2 in. (51 mm)	212241	212249	-

Casting Socks

US

Ordering information

Sold per dozen

DESCRIPTION	PRODUCT NUMBER
Casting Balloon, Large, 25 pack	199076
Casting Balloon, Medium, 25 pack	199084
Casting Balloon, Small, 10 pack	199092

Comfort Casting Socks

US

Ordering information

Sold in packages of 6

DESCRIPTION	PRODUCT NUMBER
Below Knee application	636000006



Thermoformable Composite (TFC), Sheets

SE

Product Features

- Unidirectional - Maximal strength in one direction and are good for bars, insoles and strengthening
- Bi-directional - The same strength from 0° and 90° direction and good for difficult shapes like the heel
- Sheet size: 39.4 x 23.6 in. (1000 x 600 mm)



Ordering information

COMFIL® THERMOFORMABLE COMPOSITE UNIDIRECTIONAL

DESCRIPTION	USAGE	THICKNESS	PRODUCT NUMBER
Unidirectional Carbon	Bars and other rigid parts	0.08 in. (2.0 mm)	700010101
Unidirectional Carbon & Glassfiber	Bars and rigid insole	0.1 in. (2.55 mm)	700010103
Unidirectional Carbon	Reinforcements	0.03 in. (0.7 mm)	700010126

COMFIL® THERMOFORMABLE COMPOSITE BIDIRECTIONAL

DESCRIPTION	USAGE	THICKNESS	PRODUCT NUMBER
Bi-directional Glassfiber	Opening flaps	0.03 in. (0.9 mm)	700010121
Bi-directional Glassfiber	Flexible toe off parts	0.06 in. (1.5 mm)	700010122
Bi-directional Glass-and Carbonfiber	Rigid foot plates	0.06 in. (1.5 mm)	700010123
Bi-directional Carbonfiber	Difficult shapes and strengthening	0.03 in. (0.7 mm)	700010124
Bi-directional Glass-and Carbonfiber	Rigid foot plates	0.08 in. (2.1 mm)	700010125
Bi-directional Carbonfiber	Difficult shapes and strengthening	0.02 in. (0.5 mm)	700010127
Bi-directional Glassfiber, clear	Opening flaps	0.03 in. (0.7 mm)	700010128
Bi-directional LCP	Increased flexibility	0.02 in. (0.6 mm)	700010135

For information about the Thermoformable Composite material please contact Fillauer

Thermoformable Composite (TFC), Insole

SE

Product Features

- Flexible insoles are 0.06 in. (1.5 mm) thick and are better for shaping double curves, for example around the heel
- Rigid insole is 0.1 in. (2.5 mm) thick and made to be used as a rolling insole for diabetic patients
- Made for shaping single curves and is harder to shape around the heel
- Sold in pairs



Ordering information

COMFIL® INSOLE FLEXIBLE

	EUROPE	US	PRODUCT NUMBER
Size	35 - 37	3.5 - 5	700030121
Size	38 - 40	6 - 7.5	700030122
Size	41 - 43	8 - 9	700030123
Size	44 - 46	10.5 - 12	700030124

COMFIL® INSOLE RIGID

	EUROPE	US	PRODUCT NUMBER
Size	35 - 37	3.5 - 5	700030141
Size	38 - 40	6 - 7.5	700030142
Size	41 - 43	8 - 9	700030143
Size	44 - 46	10.5 - 12	700030144




Bars

SE

Product Features

- We recommend 700020321 and 700020331 for AFO reinforcement
- We recommend 700040420 for matrix seating reinforcement
- Size 23.6 x 0.8 in. (600 x 20 mm)



Ordering information

DESCRIPTION	THICKNESS	PRODUCT NUMBER
Unidirectional, Carbon & Glass Fiber	0.1 in. (2.55 mm)	700020321
Unidirectional, Carbon Fiber	0.06 in. (1.6 mm)	700020331
Unidirectional, Carbon Fiber	0.12 in. (3.2 mm)	700040420

TFC Upright Bars & Silicone Forming Sleeve

US

Product Features

- Stronger and lighter than conventional metal uprights
- Thermoformable contouring using heat gun or oven
- Accommodates smaller radii bends
- Increased stiffness and flexibility of carbon composite
- Will not corrode unlike conventional metal uprights
- Uses innovative post pre-preg carbon material
- Optimized attachment plates for conventional metal joints
- Half oval



Ordering information

DESCRIPTION	SIZE	PRODUCT NUMBER
COMFIL® Composite Upright	3/4 x 18.5 in. (19 x 470 mm)	700600
Silicone Forming Sleeve	3/4 x 18.5 in. (19 x 470 mm)	700610
COMFIL® Composite Upright	5/8 x 18.2 in. (16 x 460 mm)	700601
Silicone Forming Sleeve	5/8 x 18.2 in. (16 x 460 mm)	700611
Flat Plate For TFC Upright 3/4 in. Bar		700620
Flat Plate For TFC Upright 5/8 in. Bar		700621

Cotton Stockinette

SE

Product Features

- Good for air evacuation for vacuum
- 3.2 in wide, 22 yd. (20 m) per role



Ordering information


DESCRIPTION	PRODUCT NUMBER
Cotton Stockinette	636508020

Teflon Sheet

SE

Product Features

- Cover in the oven when heating material up
- Dimensions 40 x 40 x 0.01 in. (1000 x 1000 x 0.25 mm)



Ordering information

DESCRIPTION	PRODUCT NUMBER
Teflon Sheet	700000010

Silicone Sheet

US



Ordering information

DESCRIPTION	USAGE	THICKNESS	SIZE	PRODUCT NUMBER
Silicone Sheet, Red	Shiny surface	0.02 in. (0.5 mm)	39.4x39.4 in. (1000x1000 mm)	700000011
Silicone Sheet, Transparent, 30 Shore A	Matt surface	0.04 in. (1.0 mm)	39.4x39.4 in. (1000x1000 mm)	700000022
Silicone Sheet, Transparent, High Gloss	Glossy surface	0.03 in. (0.8 mm)	36.2x36.2 in. (920x920mm)	600320

Vacuum Bag material

SE

Product Features

- Roll of 23 ft (10 m)
- Material: Nylon
- Thickness: 0,05 mm
- Width: 457 mm
- Maximum temperature: 246°C**



**:

 The maximum use temperature is dependent upon duration at maximum temperature and is process specific, we recommends testing prior to use.

Ordering information

DESCRIPTION	SIZE	PRODUCT NUMBER
Vacuum Bag Material, Tube	18 x 0.002 in. (457 x 0.05 mm)	700000012

Tape

SE

Product Features

- High temp. tape for use with TFC material



700000015



700000014

Ordering information

DESCRIPTION	USAGE	SIZE	PRODUCT NUMBER
Sealing Tape for Vacuum Bag	Sealing the vacuum bag	0.12 x 0.5in. x 24.9 ft. (3 x 12 mm x 7.6 m)	700000014
Tape	Fasten the release film to the plaster	1 in. x 108 ft (25 mm x 33 m)	700000015

Silicone Glove

SE

Product Features

- Sold individually
- Colors may vary

Ordering information

DESCRIPTION	PRODUCT NUMBER
TFC Silicone Glove	700070010



Release Film - sheet

SE

Product Features

- Prevent moisture from reaching the carbon fiber
- Size 47.2 x 39.4 in. (1200 x 1000 mm)

Ordering information

DESCRIPTION	THICKNESS	PRODUCT NUMBER
Release Film	0.001 in. (0.025 mm)	700000023



Silicone Bag

SE

Product Features

- For shiny surface

Ordering information

DESCRIPTION	PRODUCT NUMBER
Silicone Bag, small	700000105
Silicone Bag, medium	700000107
Silicone Bag, large	700000109





Index

Accessory Holder for Bench Carver	17	PDQ BT-4 Infrared Oven	32
Accessory Holder for Floor Carver	14	PDQ PO-2 Infrared Oven	33
Air Attachment for Floor Carver	14	PDQ Rollabout Oven	31
Alginate Impression Powder	52	Pin Wrench	43
Bars	60	Pipe Holding Fixtures	37
Bench Carver Floor Stand	17	Pneumatic Sand Drum	20
Bending Irons	40	PolyCar-C Sheets	53
Cable Swaging Tool	38	Polyester Stockinette - Lamination Tube	55
Cast Cutter	44	Polypropylene	50
Casting Socks	57	Prosthetic Bench Vise	37
Cast Mandrels	45	PVA Bags	52
Catalog symbol legend	6	Release Film - sheet	63
Cenpol	49	Rounded Blade Spatula	47
Cleani-Stump Antiseptic Wipe	51	Router Adapters	25
Comfort Casting Socks	57	Rubber Sand Cones	21
Contouring Fork and Bench Receptacle	39	Scarpas Knife	43
Contouring Irons	40	Scissors	41
Cotton Stockinette	56	Shaft Foam Arbors	19
Cutting Tool for 30 mm Tubes	44	Silicone	53
Dacron Tape	57	Silicone Bag	63
Deburring Tool	44	Silicone Glove	63
Drape Forming Adapters	29	Silicone Sheet	62
Drill Attachments 1/2"-13 Thread	18	Soft Foam Shaper	19
Dual Purpose Expansion Arbors	24	Spring Clamps	41
Dust Collection Hose Kit for 2100L	13	Surform™ Files	45
Dust Collector for Bench Carver	17	Surgical Scalpel	43
Dust Collector for Floor Carver	13	Surgical Tubing	35
Elastic Stockinette Tube, 7% Lycra	55	Tape	62
Elastic Stockinette Tube, 12,5% Lycra	56	Tapered Felt Arbor/ Core	22
Expansion Arbors	23	Teflon Sheet	61
Grinding Rasps	25	TFC Upright Bars & Silicone Forming Sleeve	60
Hex Keys	42	Thermoflex	51
Hex Screwdriver	42	Thermo Foam 4E/ Firm Foam	50
Holding Frames	29	Thermoformable Composite (TFC), Insole	59
Joint Alignment Tool	39	Thermoformable Composite (TFC), Sheets	58
Laminating Fixture LJ-100	35	Torque Screwdriver	42
Lift Assist Kit for Floor Carver	15	Torque Wrench	42
Measuring Tapes	45	Transfemoral Length Gauge	47
Original AP-ML Gauge	46	Transforming Silicone Putty	54

Transtibial Length Gauge	46
Trautman Coneiak Polishing Arbors	26
Trautman Double Drum Sander	9
Trautman Single Speed Bench Carver	16
Trautman Single Speed Floor Carver	11
Trautman Tech Station	12
Trautman Variable Speed Bench Carver	16
Trautman Variable Speed Floor Carver	10
Tycro Polishing Arbor	26
Universal Alignment Fixture	38
Universal AP-ML Gauge	46
Urethane and Felt Arbors	20
Vacuum Bag material	62
Vertical Fabricating Fixture, VFJ-100	36
Yates Clamp	41



