



Natural Definition Glove

Fillauer®

Natural Definition Glove

With the Natural Definition Glove, we take the prosthetic cosmetic glove to the next level. Using a dual layer PVC with added fibers and airbrush technique, the glove shows a deeper definition and more natural look. The glove is treated with a Micro-coated surface, making it more durable and stain resistant compared to standard PVC gloves. Environmentally responsible and sustainable, the PVC used is free from all heavy metals and phthalates which gives you a skin friendly and fully recyclable product.

Product Features

- Airbrushed, dual layer PVC with fibers for natural definition
- Available in 16 colors
- Micro-coated for durability, easy cleaning, and low friction against clothing
- Recyclable

Warranty

This product has a 6 month warranty against manufacturing defects

Natural Definition Guide

For your convenience, we made the Natural Definition Guide to help determine which color to choose for your patient.

The Guide consists of two color scales, R and Y (Red and Yellow). Each one has 8 different colors to choose from to cover most skin tones. To help your patient choose the best color, start with the color scale that most closely resembles their natural skin tone. For example, if the patient has a pinkish hue to their skin, then start with the colors on the R scale.



**Product number
611200003**



Gloves for Child Passive Hands

Style Number	Side	Approximate Age for Use	A	B	C	D	E
10901	Left	1 – 3 years	115 mm	115 mm	45 mm	145 mm	210 mm
10902	Right		115 mm	115 mm	45 mm	145 mm	210 mm
10903	Left	3 – 6 years	140 mm	115 mm	50 mm	145 mm	230 mm
10904	Right		140 mm	115 mm	50 mm	150 mm	230 mm
10905	Left	5 – 9 years	150 mm	130 mm	55 mm	160 mm	255 mm
10906	Right		150 mm	130 mm	55 mm	165 mm	255 mm

Glove Ordering Information

To order: 1 + Style Number + Color + 4 = Glove Number

Example: Style Number 10901 & Color R07 =
Glove Number 110901R074


Hand Ordering Information

5 + Style Number + 000 = Passive Hand

5 + Style Number + 001 = Reinforced Passive Hand

5 + Style Number + 002 = Super-Reinforced Passive Hand

Example: Style Number 10901 & Reinforced Passive Hand =
Hand Number 510901001



- A** Metacarpal circumference (MCP)
- B** Wrist circumference
- C** Middle finger length (palmar)
- D** Glove end circumference
- E** Length of glove (palmar)



10901
10902




10903
10904



10905
10906

Gloves for Male Passive Hands

Style Number	Side	Remarks	A	B	C	D	E
30909	Left		171 mm	143 mm	75 mm	206 mm	330 mm
30910	Right		168 mm	146 mm	75 mm	215 mm	330 mm
30911	Left		185 mm	156 mm	83 mm	215 mm	320 mm
30912	Right		185 mm	156 mm	83 mm	210 mm	320 mm
30913	Left	Prominent veins	197 mm	167 mm	90 mm	230 mm	370 mm
30914	Right		200 mm	161 mm	85 mm	225 mm	360 mm
30915	Left		190 mm	146 mm	82 mm	235 mm	370 mm
30916	Right		200 mm	168 mm	90 mm	240 mm	450 mm
30917	Left	Long fingers	214 mm	178 mm	95 mm	224 mm	400 mm
30918	Right		225 mm	186 mm	100 mm	295 mm	400 mm
30919	Left		217 mm	176 mm	92 mm	275 mm	360 mm
30920	Right		212 mm	174 mm	90 mm	240 mm	340 mm
30921	Left		222 mm	186 mm	95 mm	282 mm	360 mm
30922	Right		225 mm	177 mm	95 mm	285 mm	370 mm



- A** Metacarpal circumference (MCP)
- B** Wrist circumference
- C** Middle finger length (palmar)
- D** Glove end circumference
- E** Length of glove (palmar)

Glove Ordering Information

To order: 1 + Style Number + Color + 4 = Glove Number

Example: Style Number 30911 & Color R07 =
Glove Number 130911R074

Hand Ordering Information

5 + Style Number + 000 = Passive Hand

5 + Style Number + 001 = Reinforced Passive Hand


5 + Style Number + 002 = Super-Reinforced Passive Hand

Example: Style Number 30911 & Reinforced Passive Hand =
Hand Number 530911001



Gloves for Female Passive Hands

Style Number	Side	Remarks	A	B	C	D	E
22907	Left	Medium nails	154 mm	122 mm	70 mm	200 mm	410 mm
22908	Right		151 mm	123 mm	65 mm	200 mm	400 mm
22909	Left	Medium nails	167 mm	128 mm	72 mm	214 mm	440 mm
22910	Right		165 mm	133 mm	70 mm	215 mm	440 mm
22911	Left	Long nails	170 mm	148 mm	75 mm	205 mm	490 mm
22912	Right		170 mm	148 mm	75 mm	217 mm	490 mm
22703	Left	Long nails	174 mm	140 mm	78 mm	222 mm	470 mm
22704	Right		172 mm	138 mm	78 mm	222 mm	470 mm
22917	Left	Long nails and prominent veins	180 mm	157 mm	82 mm	245 mm	430 mm
22918	Right		184 mm	157 mm	82 mm	245 mm	430 mm
20909	Left	Unisex Short nails	162 mm	134 mm	75 mm	195 mm	320 mm
20910	Right		164 mm	142 mm	73 mm	204 mm	320 mm
20911	Left	Unisex Short Nails	170 mm	148 mm	75 mm	210 mm	380 mm
20912	Right		170 mm	148 mm	75 mm	235 mm	380 mm
20913	Left	Long fingers and nails	172 mm	146 mm	90 mm	215 mm	370 mm
20914	Right		170 mm	142 mm	90 mm	223 mm	370 mm
20915	Left	Short nails	175 mm	150 mm	84 mm	220 mm	360 mm
20916	Right		178 mm	150 mm	82 mm	199 mm	380 mm
20917	Left	Prominent veins	190 mm	162 mm	90 mm	230 mm	385 mm
20918	Right		198 mm	158 mm	85 mm	200 mm	385 mm



- A** Metacarpal circumference (MCP)
- B** Wrist circumference
- C** Middle finger length (palmar)
- D** Glove end circumference
- E** Length of glove (palmar)

Glove Ordering Information

To order: 1 + Style Number + Color + 4 = Glove Number

Example: Style Number 22911 & Color R07 =

Glove Number 122911R074

Hand Ordering Information


5 + Style Number + 000 = Passive Hand

5 + Style Number + 001 = Reinforced Passive Hand

5 + Style Number + 002 = Super-Reinforced Passive Hand

Example: Style Number 22911 & Reinforced Passive Hand =

Hand Number 522911001



Gloves for Motion Control & Otto Bock Hands

Style Number	Side	Glove	Hand	A	B	C	D	E
30713 30714	Left Right	Male Glove for 7 ¼ in. Hand	7 ¼ in.	182 mm	178 mm	63 mm	245 mm	475 mm
30715 30716	Left Right	Male Glove for 7 ¾ in. Hand	7 ¾ in.	200 mm	182 mm	68 mm	255 mm	390 mm
30719 30720	Left Right	Male Glove for 8 ¼ in. Hand	8 ¼ in.	210 mm	197 mm	73 mm	246 mm	360 mm
20713 20714	Left Right	Female Glove for 7 ¼ in. Hand	7 ¼ in.	185 mm	160 mm	75 mm	230 mm	500 mm
20715 20716	Left Right	Female Glove for 7 ¾ in. Hand	7 ¾ in.	200 mm	165 mm	75 mm	230 mm	490 mm


Gloves for Electric Hand 2000 Otto Bock

Style Number	Side	Glove	Hand	A	B	C	D	E
10711 10712	Left Right	Child Glove for 5 in. Hand	5 in.	136 mm	125 mm	41 mm	200 mm	230 mm
10713 10714	Left Right	Child Glove for 5 ½ in. Hand	5 ½ in.	147 mm	135 mm	45 mm	200 mm	230 mm
10715 10716	Left Right	Child Glove for 6 in. Hand	6 in.	162 mm	143 mm	56 mm	210 mm	250 mm
10717 10718	Left Right	Child Glove for 6 ½ in. Hand	6 ½ in.	172 mm	160 mm	64 mm	220 mm	260 mm

Glove Ordering Information

To order: 1 + Style Number + Color + 4 = Glove Number

Example: Glove number 10713 & color R07 = Glove Number 110713R074



- A** Metacarpal circumference (MCP)
- B** Wrist circumference
- C** Middle finger length (palmar)
- D** Glove end circumference
- E** Length of glove (palmar)

Fillauer®



www.fillauer.com

Fillauer LLC

2710 Amnicola Highway
Chattanooga, TN 37406
423.624.0946
Fax 423.629.7936

© 2016 Fillauer LLC
PR330/11-14-16